



Hands-Free Gross Imaging



- ◆ Mount within grossing station, over bench top, or above autopsy table
- ◆ Auto focus/Auto exposure on live preview
- ◆ Foot pedal control: one-tap field of view & camera settings
- ◆ Auto-label with case information & scale bar
- ◆ Easy touchscreen or mouse annotation & markup of image
- ◆ Quick, calibrated measurements on image
- ◆ PAXcamCONNECT™ live streaming for remote consultation
- ◆ Includes all *PAX-it!* tools for image management, annotation, processing, and sharing
- ◆ Capture directly to your case management software, to MS-Office® applications, to *PAX-it!* imaging software, or to file
- ◆ Compatible with PAXcam microscope imaging solutions for an integrated lab imaging protocol



PAXcamHD Gross Imaging System Specifications and Benefits

- ◆ **Work inside the grossing station, over a grossing bench, or over an autopsy table** without moving the specimen to a separate camera stand
- ◆ **Compact form factor** in a hermetically sealed, smooth camera body that may be wiped down
- ◆ **HD1080p live preview** of the gross image area displayed on a standard computer monitor
- ◆ **30x zoom range** driven by quick and silent motors, yielding a wide range of fields of view (FOV) at a given working distance, for example:
 - ◆ 45" wide to less than 2" wide fields of view at 38" working distance
 - ◆ 27" wide to less than 1" wide fields of view at 22" working distance
- ◆ **Auto focus** at any zoom level, with manual focus override through computer controls – *no need to touch the camera*
- ◆ **Auto exposure** with exposure compensation settings and manual exposure override through computer controls – *no need to touch the camera*
- ◆ **30 user-defined, one-tap presets** to easily store and recall any combination of the following:
 - ◆ Specific zoom, exposure and contrast levels to accommodate different specimens' size, shape, and reflectivity
 - ◆ User-designed annotation and scale bar templates for instant application to captured images
 - ◆ Defined methods for saving images from the PAXcamHD camera, for example:
 - To a currently opened case folder in *PAX-it!*
 - To a *PAX-it!* folder based on Accession Number or other user-designated sorting method
 - To a Windows file, PowerPoint® presentation, or Outlook® email, in addition to a defined *PAX-it!* location
- ◆ **Foot pedal control** for switching between the stored camera / zoom presets and for image capture – *hands-free operation!*
- ◆ **Versatile, user-friendly software interface** for use with **Touch Screen** or traditional mouse controls
- ◆ **Capture directly to *PAX-it!* Software** and other convenient locations, including: MS-Office® applications (PowerPoint®, Outlook®, Word®), to file, your hospital application, and more!
- ◆ **Capture video clips** of gross specimen handling, including time-lapse sequences
- ◆ **Includes PAXcamCONNECT™ Live Streaming** to share PAXcamHD's live preview with multiple remote locations via LAN or the internet, for interactive telepathology and consultation regarding gross specimen prep and features
- ◆ **Includes *PAX-it!* software** to seamlessly capture images into a user-defined database archive, with **touch-friendly** markup tools, reporting tools, and more:
 - ◆ Draw lines, boxes, and circles on your image manually or with easy user-defined presets
 - ◆ Add pre-formatted textual phrases to your image with one touch
 - ◆ Add a calibrated scale bar or multiple scale bars to the image
 - ◆ Draw calibrated measurement lines with immediate readout in your choice of measurement units
 - ◆ Simple tap control to drop sequential labels (A,B,C... or 1,2,3...) for serial sections or cassette designations
 - ◆ Additional annotation and measurement tools
 - ◆ Image processing tools including HDR imaging, fusion of focal planes, stitching
 - ◆ Expandable software modules, including: image analysis, network site licensing, motorized stage control, and more
- ◆ **Supplemental ring light illumination** is available for setups where ambient light is insufficient
- ◆ **Optional camera mounts** are available, including flexible arms for mounting within a grossing station, and fixed or articulating arms for bench top installations

(Microsoft, Word, Windows, Office SBE, Outlook, Excel and PowerPoint are registered trademarks of Microsoft Corporation)

MIS MIS Worldwide Headquarters
707 N. Iowa Avenue
Villa Park, IL 60181 USA
+1.630.279.4000 tel. +1.630.279.4050 fax.

PAXcam™
www.paxcam.com

PAX-it!
www.paxit.com