

HISTO-LOGIC[®]

No reader should utilize materials and/or undertake procedures discussed in HISTO-LOGIC articles unless the reader, by reason of education, training and experience, has a complete understanding of the chemical and physical properties of all materials to be utilized and the function of each material by which any procedure is accomplished.

Editor, Lee G. Luna, D. Litt., H.T. (ASCP)

Technical Bulletin for Histotechnology
Published: January, April, July, October

Volume XII, No. 2 - April, 1982

In Defense of Cover Glasses

Gary Gill, CT (ASCP), CMIAC, and

John Frost, M.D., FCAP, FIAC

Division of Cytopathology

The Johns Hopkins University School of Medicine

Baltimore, Maryland 21205

The recent recommendation for the use of liquid substitutes for cover glasses (Buxsa M: Liquid coverslipping media. *Histo-Logic* 11(4):167-168, 1981) compels me to present several observations about cover glasses that argue against the use of liquid substitutes. First made available commercially in 1840 by Chance Brothers of Birmingham, England,¹ cover slips should be regarded as more than just a "pane in the glass," as it were, that is used in a costly, time-consuming procedure and that drives people to look for a less costly and less tedious alternative. After all, there must be some valid reason why, ounce for ounce, cover glasses cost about 15 times more than micro slides.

Cover glasses are manufactured to meet well-defined optical and physical specifications (e.g., refractive index of 1.523 ± 0.005 , dispersion value of 52.0, thickness, plane parallelism) that are related to the design of microscope objectives and that visibly influence image quality.^{2,3} Indeed, cover glasses alter the path of light so substantially that lens designers consider them as the first lens in all microscope objectives! By extension, coverslipping can be regarded as the last step in an objective's assembly. Seen in this light, coverslipping deserves the best materials and methods.

Given a well-prepared cell spread or tissue section and a clean microscope adjusted to provide Kohler illumination, the correct thickness of cover glass and mounting medium will reveal crisp clear images. Too thick or too thin cover glass and mounting medium layers introduce spherical aberration that overall lowers image contrast, thereby producing images described variously as blurry, hazy, milky, watery, soft focus, out-of-focus, unfocusable, cloudy, flat, foggy, fuzzy, and washed-out. Edges and boundaries of cellular details become unsharp. Images are degraded in direct proportion to the degree of deviation of the actual thickness of cover glass and mounting medium from the intended thickness. The degree of image degradation can range from barely perceptible to ruinous.

Microscope objectives tolerate deviations from the design-specified thickness of cover glass and mounting medium as a function of their numerical aperture (NA). Numerical apertures range from about 0.04 to 0.95 for dry objectives, 1.0 to 1.3 for oil immersion ones, and are a measure of the light-gathering capacity of an objective — the higher the NA of an objective, the greater its resolving power. 10X objectives have NAs of 0.25 to 0.32 and are insensitive to cover glass thickness. Indeed, a micro preparation can be examined upside-down, with the slide serving as a thick cover glass, without noticeably degrading the image. On the other hand, 0.65 NA achromat objectives (i.e., those least corrected for optical aberrations) tolerate deviations only within $\pm 15 \mu\text{m}$ of 0.18 mm (American- and English-made microscopes) or 0.17 mm

(European- and Japanese-made microscopes), and 0.95 NA apochromat objectives (i.e., those most highly corrected for optical aberrations) tolerate deviations only to within $\pm 3 \mu\text{m}$! Indeed, the latter objectives would be impractical to use were it not for their having an adjustable correction collar that alters the spacing among the lens elements to compensate for cover glass and mounting medium thicknesses of 0.11 to 0.23 mm.

The thickness of cover glass recommended most often is the No. 1½. No. 1½ cover glasses range in thickness from 0.16 or 0.17 to 0.19 mm.^{2,3} This range brackets the 0.17 and 0.18 engraved on some microscope objectives, numbers that are usually interpreted as indicating the thickness of cover glass only, no mounting medium, for which the objective should be used. No. 1½ cover glasses are correct only when there is little or no mounting medium between the cells or tissue and the cover glass. Examples include: (1) monolayer cell spreads that are heated on a hot plate (i.e., "cooked") to quickly evaporate the solvent of the mounting medium, thus leaving behind only a very thin layer of hardened mounting medium; (2) cells spread directly onto a cover glass (e.g., a blood film); (3) cells grown in culture on a cover glass; and (4) cells collected on a Nuclepore filter that is to be dissolved cell-side-down on a cover glass. These four examples are uncommon.

In routine practice, the No. 1 thickness cover glass is more correct. No. 1 cover glasses range from 0.13 to 0.17 mm in thickness. While thinner than usually recommended, No. 1 cover glasses allow for mounting medium, the consequences of which are curiously overlooked in most recommendations (even though mounting medium interacts with light as though it were glass). The thickness of mounting medium can often be substantial, exceeding by several times the most relaxed tolerance limits⁴ (i.e., $\pm 15 \mu\text{m}$ for 0.65 NA achromats). Therefore, it is wise to use No. 1 thickness cover glasses and to apply as little mounting medium as will fully cover the specimen without subsequent retraction under the edge of the cover glass. The combined thickness of a No. 1 cover glass and mounting medium will more closely and more often equal the specified thicknesses of 0.17 or 0.18 mm "cover glass" for which objectives are corrected than a No. 1½ cover glass and mounting medium ever can.

The ruinous spherical aberration introduced by large deviations in cover glass and mounting medium thickness from the specified thickness, for 40X objectives especially, can be offset substantially by closing the aperture diaphragm of the substage condenser more than usual when establishing Kohler illumination in a microscope. The aperture diaphragm is ordinarily focused in the back focal plane of the objective and so determines the objective's working numerical aperture. Closing the aperture diaphragm, though never completely, reduces the objective's numerical aperture and lessens its sensitivity to cover glass thickness. In practice, simply close the aperture diaphragm until the best contrast is produced.

The foregoing is but a minimum introduction to some important aspects of cover glasses: a thorough discussion would require many more pages. Undoubtedly, not all

readers will have been convinced to continue to use cover glasses instead of liquid substitutes. To them in particular, and to all in general, it must be emphasized that the beneficial effects of cover glasses, as well as the detrimental effects of liquid substitutes, can be masked by poorly preserved specimens, pale staining, dirty microscopes, improperly adjusted microscopic illumination (e.g., the low contrast produced by spherical aberration is visually indistinguishable from that produced by glare when an aperture diaphragm is opened excessively or dirty lenses scatter light), or are observed only under a 10X objective. Under professionally controlled circumstances, the difference in the quality of images made possible by cover glasses and liquid substitutes is real, visible, and substantial. If good image quality is among one's goals for specimen preparation, then only properly selected and applied cover glasses can make it possible.

About the best that can be said for liquid substitutes for cover glasses is that they are inexpensive, fast-drying glues that are quick to use. In contrast to cover glasses and most mounting media, liquid substitutes: (1) produce greater volumes of potentially toxic vapors; (2) do not produce a coating of uniform thickness; (3) produce a wavy layer; (4) have uncontrollable thicknesses, often being thinner than is best, that will produce spherical aberration and hazy images; (5) have a low refractive index relative to that of cells and so produce a refractile appearance in cells; (6) can leave thick areas of cells uncovered; (7) may increase stain fading (the inclusion of an anti-oxidant does not ensure the stability of all dyes); (8) scratch easily; and (9) are not readily cleaned of immersion oil by xylene.

"If liquid substitutes for cover glasses are so unsuitable," it might be asked, "why have they been introduced to the microtechniques market?" Although only the manufacturers of these substitutes can answer this question, I suspect that someone saw a business opportunity in providing a product that would quickly and inexpensively cover a micro preparation with a thin protective film of rapid drying transparent medium. Without knowing the information in the cited references and others, a manufacturer might never stop to ask whether a cover glass has specific properties that a liquid substitute should have as well. Whatever their evolutionary path from concept to market, liquid coverslipping media can not truly substitute for cover glasses; the laws of optics will not allow it. Using liquid substitutes for cover glasses means accepting limitations in the quality of one's work. So while much routine microscopic work can be done satisfactorily with less than optimal micro preparations, we should nonetheless strive to maintain the high standards of performance that are the hallmark of professionals and use cover glasses.

All things considered, liquid substitutes offer advantages only to the persons who must coverslip and to any administrator who does not have responsibility for quality of patient care or malpractice suits arising from less-than-accurate microscopic diagnoses. Conventional cover glasses benefit microscopists who must see highly resolved detail, especially pathologists with patients' health and lives at stake, by allowing them to see cells and tissues more clearly, and thereby increase the likelihood of their interpreting pathologic changes accurately. In the total picture, it is the needs of the microscopist and the patient that must be served. Liquid substitutes will not do, for there is more to a cover glass than meets the eye and the pocketbook.

References:

1. Braccagirdle, B.: *The History of Microtechnique*. Cornell University Press, Ithaca, 1978, p. 113.
2. ASTM Committee E-25 on Microscopy: Standard Specification for Cover Glasses and Glass Slides for Use in Microscopy - Designation E211-70, in *Annual Book of ASTM Standards*, American Society for Testing and Materials, Philadelphia, 1970.

3. General Services Administration: Interim Federal Specification - Cover Glass, Slide, NNN-C-001434A (GSA-FSS), GSA, Washington, D.C., Jan. 8, 1971.
4. Settingington, R.: The specification of a standard microscope cover glass. *J. Roy. Micro. Soc.*, 73:69-76, 1953.
5. Aumonier, F. J., and Settingington, R.: Some Notes on the Mounting of Histological Sections. *Proc. Roy. Micro. Soc.*, 2:428-429, 1967.



Bis-Chloro-Methyl Ether (BCME): A Response

Gary Gill, CT (ASCP), CMIAC
Department of Cytopathology
The Johns Hopkins Hospital
Baltimore, Maryland 21205

Janet Maass' concern regarding the possible formation of bis-chloro-methyl ether (BCME), a potent carcinogen, as a result of formaldehyde being mixed with hydrochloric acid is understandable (*Histo-Logic* 10(4):152, 1980). Fortunately, however, published data indicate that BCME formation is highly unlikely under the conditions of use ordinarily encountered in safely operated histology laboratories.¹

Whether BCME will result from exposure of formaldehyde to HCl depends on the concentrations of each in parts per million (ppm) in air, not in solution. When the concentrations of both chemicals are kept below their respective threshold limit values - time weighted average (TLVs-TWA), no detectable BCME is formed.¹ (TLV-TWA is "the time weighted average concentration for a normal 8-hour workday or 40-hour workweek, to which nearly all workers may be repeatedly exposed, day after day, without adverse effect.")² The TLV-TWA for formaldehyde is 2 ppm; for hydrochloric acid (hydrogen chloride), 5 ppm. The TLV-TWA for BCME is 0.001 ppm or 1 part per billion (ppb). Such low values for all three chemicals place them in the highly toxic category. For comparison, the TLVs-TWA for xylene and ethyl alcohol are 100 ppm (upper limit of the highly toxic range) and 1,000 ppm (slightly toxic), respectively.³

Less than 0.5 ppb of BCME is formed from 20 ppm each of formaldehyde and HCl under a variety of conditions of humidity, temperature and UV radiation.¹ In other words, even when the gaseous concentrations of formaldehyde and HCl exceed their TLV-TWA by factors of 10 and 4, respectively, the amount of BCME produced within 12 to 24 hours is within safe limits.¹ Concentrations of formaldehyde and HCl this high are unlikely to occur, for the noxious vapors would alert the laboratory worker to take corrective action long before reaching such olfactorily intolerable levels.

While some data conclusively demonstrate that BCME is formed from gaseous formaldehyde and HCl, other data suggest that acidic solution (pH unspecified) which contains chloride ion and formaldehyde should be considered a potential source of BCME.⁴ Whether B-5 fixative (sodium acetate-mercuric chloride-formalin: 9 parts distilled water, 1 part 40% formaldehyde, 6% HgCl₂, and 1.25% NaC₂H₃O₂ [anhydrous], pH 5.8 to 6)⁵ is sufficiently acidic to produce BCME is uncertain, though it seems highly unlikely.

Knowing the TLVs-TWA of formaldehyde, hydrogen chloride and other toxic vapors are useless when we do not know the actual concentrations of vapors to which we are exposed. In the absence of such information, a useful rule-of-thumb is that the concentration of any odorous toxic chemical is excessive when we can smell it. The minimum concentration at which the odor quality of a compound under study can be described is its recognition odor threshold concentration (ROTC). A chemical's ROTC is usually substantially

National Society for Histotechnology Symposium/Convention

Boston, Massachusetts

September 13-17, 1982

The Eighth Annual Symposium/Convention of the National Society for Histotechnology will be conducted at the Sheraton Hotel in Boston. The enclosed program is complete with hotel reservation card and registration form.

HOTEL RESERVATIONS: We request you make hotel accommodations early, as all rooms blocked for NSH will be released one month prior to meeting date. **MAIL HOTEL RESERVATION DIRECTLY TO:** Sheraton, Prudential Center, Boston, Massachusetts 02199; 617/236-2000.

SYMPOSIUM REGISTRATION: Registration form may be photocopied if more than one individual from the same activity desires to attend. To avoid delays and unnecessary complications, registrations **AWAITING FUND APPROVAL** will be accepted and held in abeyance until final commitment is received. Please include a note to this effect on your registration form. For workshops stipulating a **LIMITED ATTENDANCE**, please indicate 1st and 2nd choice. Your first choice will be honored as long as openings remain. Workshop quotas will be adhered to. **LATE REGISTRATION CHARGE** of \$10 is applicable to registrations received after September 7th. Attendees registering upon arrival to the meeting will be assessed the \$10 late fee. **CANADIAN AND FOREIGN** attendees should remit fees in U.S. currency. Checks and money orders are payable to: National Society for Histotechnology. **MAIL SYMPOSIUM REGISTRATION AND FEES DIRECTLY TO:** NSH, P.O. Box 36, LANHAM, MARYLAND 20706. For registration assistance, please call Roberta Mosedale, NSH Office, 301/552-9678.

REIMBURSEMENT POLICY: NSH will reimburse registration fees upon receipt of cancellation prior to September 7th. Fees not refundable after this date. After arrival to the meeting, **NO REFUNDS** will be made for unattended workshops, sessions or banquet. Banquet tickets may be purchased at the registration desk, but are not refundable. Monetary difference cannot be refunded when changing workshop attendance after arrival to the meeting.

NSH/Thomas Edison Program Schedule

COURSE REVIEWS: Students interested in review sessions should register for individual scheduled workshops below. Workshop registration is free to students formally enrolled in Thomas Edison College. **STUDENT MUST SEND PROOF OF ENROLLMENT WITH REGISTRATION FORM.**

HUMAN MICROSCOPIC ANATOMY

(Freida Carson, Ph.D., & John Ryan, HTL [ASCP])
MONDAY, WORKSHOP NO. 2 - 8:30 AM - 4:30 PM

Human Microscopic Anatomy workshop/review will be continued Tuesday morning for those planning to take this examination during the convention week. Review will allow informal discussion with faculty member and an opportunity for examinee to identify and review weak areas. Examinees are encouraged to attend both sessions.

INTRODUCTORY HISTOTECHNOLOGY/HISTOCHEMISTRY

(Richard Schroeder, M.A. & Jules Elias, M.A.)
TUESDAY, WORKSHOP NO. 13 - 8:30 AM - 4:30 PM

EXAMINATION SCHEDULE: WEDNESDAY, THURSDAY & FRIDAY - 7-9 AM

NOTE: Registrants planning to take Thomas Edison examination(s) during convention week must *pre-register* by sending NSH your name and address with examination title(s) you will be taking. This information may be included with your convention registration or submitted separately, but is more important for administrative purposes. Examination fees may be paid at time of administration; however, advanced payment is preferred and should be mailed to: NSH, P.O. Box 36, Lanham, Maryland 20706.

National Society For Histotechnology Symposium/Convention

September 13-17, 1982

Hotel Reservation Form

SHERATON-BOSTON HOTEL & TOWERS

MASS ROOM TAX 5.7%

To assist us in accommodating your guest room request please follow these guidelines:

A) Group reservation requests will be accepted until four weeks prior to arrival. Thereafter, they will be handled on a space available basis.

B) Every effort will be made to either honor the rate requested or assign the nearest available rate.

C) Check out time is 1 pm; therefore, your room will not be available for occupancy until 3 pm. Please make your travel arrangements to coincide with these times.

D) Only one reservation card necessary per room.

ARRIVAL DAY/DATE _____ ARR. TIME _____ DEPARTURE DATE _____

GTD. _____ after 6:00 pm arrival must be accompanied by first night's rental or major credit card to guarantee the room for late arrival.

CREDIT CARD NAME _____

CREDIT CARD NUMBER _____

YOUR NAME _____

COMPANY _____

ADDRESS _____

CITY _____ STATE _____ ZIP _____

SHARING WITH _____

	NORTH BLDG		SOUTH BLDG		LUXURIOUS SHERATON-TOWERS
SINGLE	<input type="checkbox"/> 81	<input type="checkbox"/> 87	<input type="checkbox"/> 87		<input type="checkbox"/> 95
DOUBLE	<input type="checkbox"/> 93	<input type="checkbox"/> 100	<input type="checkbox"/> 100		<input type="checkbox"/> 106

MAIL DIRECTLY TO: SHERATON-BOSTON HOTEL, Prudential Center, Boston, MA 02199

WORKSHOPS

MONDAY, SEPTEMBER 13, 1982

No. 1 IMMUNOPEROXIDASE - THEORY, METHODS & CLINICAL APPLICATION (Elaine Kelen, HT, ASCP) LIMIT 30

This workshop will begin with an introductory lecture on basic immunology, immunohistochemistry and the various immunoperoxidase methods. Each registrant will also participate in an actual staining procedure. The session will conclude with a kiosk-style presentation of stained sections. Detailed handouts will be provided.

No. 2 HUMAN MICROSCOPIC ANATOMY (Intermediate to Advanced) (Doris Carson, Ph.D., & John Ryan, M.B.A., HT, ASCP)

Workshop will emphasize basic tissues and their organization (into the different organs). Visual identification, function and some special staining correlation will be included. Participants will be given a self-assessment test, and handout material will be in a workbook format to be filled in during the workshop. Workshop is recommended for Thomas Edison students.

There will be a special review on Tuesday morning for students taking the Thomas Edison examination. (Not to be taken in lieu of the workshop.)

No. 3 HISTOTECHS LOOK AT THE ENDOCRINE SYSTEM (Basic to Intermediate) (Elizabeth Swick, HT, ASCP & Elizabeth Cardillo, HT, ASCP)

This workshop will consist of lectures on the gross and microscopic anatomy and physiology of the organs of the Endocrine System, and stains used to demonstrate the tissue elements. Participants will perform selected special stains including one immunoperoxidase technique, and will be able to recognize properly stained slides. Glands covered in depth will be the pituitary, thyroid, pancreas and adrenals.

No. 4 PHOTOGRAPHY IN THE LABORATORY (Robert Kersch)

Workshop will include lecture and demonstrations in the areas of both photomicrography and gross specimen photography. Information will cover all phases of instrumentation, types of film, technique and trouble shooting routine problems.

- | | |
|------------------------|----------------------------|
| PHOTOMICROGRAPHY | GROSS SPECIMEN PHOTOGRAPHY |
| a. The microscope | a. The camera |
| b. The specimen | b. Lighting |
| c. The camera | c. Film and Processing |
| d. Film and processing | d. Projection |
| e. Projection | |

Slide presentation covering both areas of photography will be presented. Question and answer sessions will cover specific problems presented by attendees.

No. 5 SPECIAL STAINS FOR DEMONSTRATING MICRO-ORGANISMS (Charles Charafian, B.A., HT, ASCP) LIMIT 40

Workshop will present an overview of the various types of micro-organisms and the recommended staining methods used for their identification. Participants will have the opportunity to perform a rapid Gram's methanamine-silver nitrate (GMS) method for fungi and Pneumocystis carinii, a modified Brown-Bren method for Gram-positive and Gram-negative bacteria and Truant's fluorescent method for acid fast bacteria.

No. 6 PLASTIC SECTIONS - AN OFTEN NEEDED SUPPLEMENT TO PARAFFIN (Patricia Norman, M.A., HT, ASCP) LIMIT 25

Participants, through audio visual aids and hands-on experience, will become familiar with all phases of making plastic sections for use with the light microscope. Workshop will also include several case studies demonstrating how plastic, in conjunction with paraffin sections, is invaluable in making a complete diagnosis.

No. 7 TECHNIQUES FOR ARGENTAFFIN AND ARGYROPHIL CELLS (Barbara Hewitt, HT, ASCP) LIMIT 50

The lecture will include the identification of both argentaffin and argyrophil cells on 2x2 kodachrome slides using the following staining methods:

- Schmorl's Method for Reducing Substances
- A Modified Papanicolaou's Argyrophil Method
- Fontana-Masson Method for Argentaffin

These three methods are the stains the participants will perform in the wet workshop.

No. 8 INSTRUCTIONAL TECHNIQUES USED IN TEACHING (Basic) (Jules Elias, M.A., & Janet Moss, HT, ASCP)

This workshop is designed for talent conferences. Workshop to include small group instruction, writing goals for a curriculum, effective communication, writing behavioral objective, preparing activity lists for instruction and making simple visual aids. Recommended for individuals interested in using talent.

No. 9 HOW TO PLAN AN EXPERIMENT, WRITE A SCIENTIFIC PAPER AND PRESENT DATA AT A SCIENTIFIC MEETING (Charles Culling, FIMSLC, MRC, PATH)

This will be a hands-on workshop where participants will actually write a paper for publication. One participant will be selected and placed on the Scientific Session Program for Friday, to present their paper written during this workshop.

No. 10 MICROSCOPIC ANATOMY, PART I - BASIC (Roy Kerson, M.D., & Henry Albarado, HT, ASCP) LIMIT 100

An introduction to cell types, their variations, modifications, patterns, arrangements. Brief mention will be made of function and delineation of cell types by special stains. Method used for teaching will consist of study of a set of microscopic slides and projection of kodachromes with lecture. Slide sets and handouts will be kept by registrants.

No. 11 TECHNICAL AND DIAGNOSTIC ASPECTS OF THE LEDER'S STAIN (R. Peter Wittmann, M.D., & Barbara Tereola, HT, ASCP) LIMIT 50

The Leder's Stain is a histochemical, enzyme selective stain used for the demonstration of neutralophilic myeloid cells, tissue mast cells, and Chediak-Higashi inclusions. The primary significance of this enzyme stain is that it can be performed on formalin fixed paraffin embedded tissue. This workshop will give participants the opportunity to perform the more intricate steps in the Leder's stain. Information on the diagnostic significance of the stain will also be given.

No. 12 MUSCLE BIOPSY WORK-UP - HOW AND WHY? (Louis Sivas, B.S., HT, ASCP) LIMIT 75

An in-depth look at the muscle biopsy, including normal histology of muscle, how biopsies are taken, the role of the histotechnician in biopsy preparation, enzyme histochemistry and a discussion of musculoskeletal disease.

No. 13 NOMINAL GROUPING PROCESS APPROACH TO PROBLEM SOLVING (Joyce Moore, HT, ASCP)

Nominal grouping process allows a group to examine its own needs and interests and develop programs or solutions based on these findings. This workshop is designed to provide the student with a tool for identifying elements of a problem situation, identifying elements of the solution and establish priorities. The student will be taken through the steps in the Nominal Group Process. Most problems in a laboratory are created by the entire staff. Nominal grouping is an approach shared by all staff in creating a solution to the problems.

TUESDAY, SEPTEMBER 14, 1982

No. 14 INTRODUCTORY HISTOTECHNOLOGY HISTOCHEMISTRY (Basic) (Richard Schneider, M.A., & Jules Elias, M.A.)

Workshop is designed as an introduction and refresher for the discipline of histotechnology. Concepts of fixation, tissue preparation, sectioning and staining will be presented. Staining procedures utilized as routine in histopathology laboratories, i.e., carbohydrates, lipids, proteins, minerals, bacteria are discussed. The histochemistry program provides the participant with a more in-depth understanding of routine and sophisticated procedure mechanisms. Crystal and cryogenic techniques are discussed. Workshop is recommended for Thomas Edison students.

No. 15 THE USE OF HISTOCHEMISTRY IN DIAGNOSTIC PATHOLOGY: AN IMPORTANT ROLE OF THE HISTOTECHNOLOGIST IN PATIENT CARE (Hugh McAlister, M.D., C.M., D.A.)

The primary object of this workshop is to discuss the utilization of special histochemical stains in establishing diagnoses. This role of the histotechnologist in patient care will be illustrated with actual case presentations, as appropriate. Emphasis will be placed on the practical approach to diagnosis of various disease categories.

No. 16 THE HISTOLOGY AND HISTOTECHNOLOGY OF THE ENDOCRINE PANCREAS (Robert Scully, B.S., HT, ASCP & Carl Embick, HT, ASCP) LIMIT 35

Workshop is aimed to present current knowledge and concepts of the content and structure of the Endocrine Pancreas, (Islets of Langerhans), with particular reference to the changes that occur in both Zollinger-Ellison and Verner-Morrison syndromes. Patient methodologies will be presented and photographs and histological sections will be reviewed. Participants will be able to do their own immunoperoxidase.

No. 17 SEVEN KEYS TO SUCCESS IN MANAGEMENT FOR HEALTH CARE PROFESSIONALS (Jack Broderick & Phyllis Russell, B.A.) LIMIT 50

- Performance based on strength and reality.
- Managerial styles that work.
- Time is your best management resource.
- Negotiation from a position of strength and leadership.
- A win-win communications system.
- You as a "change agent" in the 80's.
- Your career development handwagons.

You gain insight into the necessity of developing skills and techniques for you and the people that work with you. You should improve your ability to create a work environment that really works, that increases worker response and department productivity. Job stress is reduced to manageable levels, better results are achieved, and employee motivation and job attitude are enhanced.

No. 18 MICROSCOPIC ANATOMY, PART II - ADVANCED (Roy Kerson, M.D., & Henry Albarado, HT, ASCP) LIMIT 100

A continuation of the analysis of cells. Emphasis will be placed on the arrangement of cells into tissues. Patterns of arrangement and techniques for analyzing patterns will be discussed. Unifying principles for histologic identification will be presented. Lecture will be presented with kodachromes. Slide sets will be kept by registrants.

No. 19 CRYTOMY TECHNOLOGY (Frank Aronson, B.A.)


This workshop is designed to teach all those individuals, technicians, technologists, and researchers both experienced and non-experienced who wish to become more competent in the field of cryotomy. The program contents will include a multitude of subjects related to all phases of the art. A brief history of frozen sectioning, the technique and instrumentation from its inception to present day. The main contents of the workshop will include such areas as detailed information on sectioning, orientation of specimen, cutting and section evaluation. The problems encountered the reasons for and the solutions to. Anti-roll plate technology, adjustment, problems, solutions and repairs. Cryostat knives, types, sharpening, angle adjustments. The cryostat, purchasing, cleaning and maintenance. The problems associated with the cryostat and cryostat evaluation. Staining methods for routine surgical and other related tissues. Histochemical reactions, immunofluorescence and immunoperoxidase techniques. Kidney and muscle biopsy technology will be discussed at length. Methods for handling of specimens shipping, freezing, sectioning and staining. Frozen section artifacts, mounting media, coverslipping, slide coatings and transport media will be discussed, in addition to numerous tricks of the trade. A cryostat will be available for demonstration and instruction for those who desire additional help.

No. 20 PRINCIPLES AND PRACTICE OF TISSUE PROCESSING (Paul Moffitt, B.S., HT, ASCP)

This workshop will provide a more thorough understanding of tissue processing by considering the various stages of processing at a molecular level. Principally for histotechnologists who do not have a background in basic physics and chemistry. Part I will consist of a review of basic physical/chemical principles which are involved in tissue processing, such as molecular motion, diffusion, density, viscosity, miscibility, convection, etc. Part II will consider specific features of automated tissue processors (vacuum, temperature, agitation, etc.) in terms of their effects on the physical/chemical processes previously discussed. Part III will offer participants an opportunity to examine several different models of automated tissue processing equipment, and to evaluate their respective features.

20% PRICE REDUCTION . . . \$4,150

complete system



Now histology's most popular Tissue Embedding System is available at a new, reduced suggested list price. And that's not all! With each system purchased you'll get a valuable bonus offer . . . a **full case** of Tissue-Tek® III Uni-Cassette™ Embedding Cassettes.

Price reduced over \$1,000

Overwhelmingly positive response to Tissue-Tek III TEC System (4584) has made it possible for us to reduce its suggested list price (from \$5,200). Add the system to your histology laboratory and take advantage of a price reduction exceeding \$1,000 effective May 1, 1982.

Get a free case of Tissue-Tek III Uni-Cassettes, too!

Purchase a TEC System and we'll send you **one full case** of Tissue-Tek III Uni-Cassette Process/Embedding Cassettes (4170), used in leading histology laboratories worldwide for over two years. Tissue-Tek III Uni-Cassettes provide convenient, one-time specimen identification from grossing to filing. And Tissue-Tek III Uni-Cassettes feature unique dual snap-latch and hinge-lock closures preventing premature separation of lid and base during processing.

Why invest in the Tissue-Tek III TEC System?

Performance, dependability and convenience make the TEC System the best you can buy . . . at any price. It's an investment in efficiency and reliability. The complete system includes separate Thermal Console (4585), Paraffin Dispensing Console (4586) and Cryo Console (4587).

Step-saving uni-plane design puts paraffin dispenser, heated work area, thermo-electric cold spot, large, refrigerated cold plate, even base mold storage area all within easy reach. This lets you move freely and smoothly from one console to the next.

Tissue-Tek III TEC System starts working for you the day it enters your lab

Reliable, solid state electronics contribute to consistently dependable performance, no matter how busy your laboratory gets. Day after day. And this performance continues for years to come.

Ask for Tissue-Tek III TEC System

Contact your distributor today and place your order for Tissue-Tek III Tissue Embedding Console System (4584). And when you order convenient Uni-Cassettes Process/Embedding Cassettes, specify the Tissue-Tek name. Accept no substitutes. Tissue-Tek equipment and accessories are available only through American Scientific Products, Curtin Matheson Scientific and VWR Scientific. Please mention this announcement when you place your order.

Write to us for more information

Return the coupon below and we'll be happy to send you more information about the TEC System.

Offer and price good only in U.S.A.
Offer expires 12/31/82

☐ Please send me more information about Tissue-Tek III Tec System (4584)

My supplier is ☐ American Scientific Products
☐ Curtin Matheson ☐ VWR Scientific

Name

Company/Institution

Street Address

City State Zip

Mail to:

LAB-TEK DIVISION, Miles Laboratories, Inc., 30W475 North Aurora Road, Naperville, IL 60566

FREE UNI-CASSETTE PROCESS/EMBEDDING CASSETTES!

See reverse side for
details . . .



FREE BONUS OFFER!

Purchase a complete Tissue-Tek III Tissue Embedding Console System (4584) and we'll send you a FULL CASE of our popular 4170 Uni-Cassettes. To receive your free case of Uni-Cassettes, simply mail us the Uni-Cassette Bonus Certificate packed with your new TEC System. (Offer expires 12/31/82.)



TISSUE-TEK III UNI-CASSETTE™ SYSTEM

... the best idea in cassettes yet! Check these unique advantages:

- Proven reliability — used by satisfied customers worldwide for over two years!
- Permanent, complete specimen identification from grossing to filing.
- Three large, specially textured writing surfaces.
- Available in five colors: white, blue, pink, green, yellow.
- Hinged lid snaps closed in one easy motion.
- Dual safety closure: snap latch and hinge lock.
- Lid stays attached to base until you're ready to discard it.
- Built in stability through cross-shaped base reinforcement.



4170

Uni-Cassettes

Slotted design Uni-Cassette insures thorough and rapid exchange of processing fluids. The Uni-Cassette's attached lid remains secure during processing because of a dual closure. When processing is complete, simply snap the lid off and discard it.

Uni-Cassettes are available in white (4170), supplied 1,500 per case. Also available in four colors for easy identification: Yellow (4154), Pink (4155), Green (4156), Blue (4157). Colors supplied 500 per case.

Biopsy Uni-Cassettes


Tissue-Tek III Biopsy Uni-Cassettes (4172) are similar to Uni-Cassettes, except for tiny 1mm pores instead of slots for the exchange of processing fluids. These pores eliminate the need to wrap most specimens and greatly reduce specimen loss. Distinctive gray color. Supplied 500 per case.



4172



4173



4166

Mega-Cassettes™ and Reusable Base Mold

For thicker, larger specimens, use Tissue-Tek III Mega-Cassette (4173). Large, textured writing surfaces eliminate the need for paper tags and help prevent identification errors. In white only. 30x24x11 mm. Note: Use with base mold 4166, below. Supplied 750 per case.

Metal base mold (4166) has deeper configuration to accommodate the Mega-Cassette, 31x23x13.5 mm. Supplied 6 per case.

☐ Please send me your complete Tissue-Tek III Histology Systems Catalog and more information about the Tissue-Tek III Uni-Cassette System

Name _____

Company/Institution _____

Street Address _____

City _____ State _____ Zip _____

My supplier is:

☐ American Scientific Products ☐ Curtin Matheson

☐ VWR Scientific

Mail to:

LAB-TEK DIVISION, Miles Laboratories, Inc.
30W475 North Aurora Road, Naperville, IL 60566

LAB-TEK DIVISION, Miles Laboratories, Inc.
Naperville, IL 60566

Lab-Tek
Division



WATER SOLUBLE EMBEDDING MEDIA FOR LIGHT MICROSCOPY

(Jane Woods, *Product Specialist*)

LIMIT 25

Recently there has been an increased interest in the use of plastic embedding for many histopathologic techniques. One of the most versatile and widely used plastics for histology is glycol methacrylate. This workshop will focus on the advantages and disadvantages of using water soluble plastic as compared with the more traditional embedding medium, paraffin. In the course of the workshop, the participants will experience the relative ease of embedding, sectioning and staining with plastics for light microscopy.

No. 24
CONNECTIVE TISSUES - FAMILY TREE AND STAINING CHARACTERISTICS(Lee D. Liss, *HT (ASCP)*)

Purpose of this course is to familiarize participants with the entire family tree of the connective tissue system. This will include the development from Mesoderm to Mesenchyme, primary, connective tissue proper and supporting connective tissues. Included will be the forms of connective tissue: i.e., fibroblast, fibrocyte, collagen, elastin, reticulin, osseous, ligaments, cartilage, mucous, adipose, tendons and blood which is a form of specialized connective tissue. The function of connective tissue will also be included. This will be followed by an in-depth discussion of the preferred fixatives, processing method, and staining procedures for the demonstration of all connective tissues. Course will be conducted with the aid of charts and photomicrographs.

No. 25
EXAMINATION DEVELOPMENT: CERTIFICATION CLASSROOM(Mary Lynn, *Ph.D.*)

LIMIT 50

The procedures for writing and evaluating multiple choice items for certification (BOH) and classroom use is the focus of this workshop. The following topics will be covered:

(1) Writing item stems and responses; (2) evaluating item stems and responses; (3) taxonomy levels in developing examination items; (4) statistical analysis of examinations; (5) examinee reporting mechanisms including score profiles and item descriptions. Participants should plan to bring five multiple choice test questions they have developed and a histology reference book. Participants will work together to develop and critique examination items as well as learn basic theories of examination construction utilized by the Board of Registry.

No. 26
DESIGN FOR HEALTH(Norman V. Steere, *P.E.*, & Diane Barrie, *HTL (ASCP)*)

This half-day workshop will address the issue of laboratory design as it relates to the Histology Laboratory. A special area of concern is ventilation, and we will look at methods of controlling fumes which are inherent in Histology, i.e., formaldehyde and xylene.

No. 27
PROCEEDINGS ON THE FOURTH BASIC SCIENCE WORKSHOP IN HISTOLOGY(A. R. Villanueva, *Ph.D.*, & J. Luis Ellis, *M.A.*)

Workshop will include the following lecture topics:

- Bone Marrow Cells - Identification, Morphology and Pathology (Rocco Letera, *M.D.*)
- Progress in Staining Decalcified Bone and Teeth (John Hess, *D.D.S.*)
- Immunohistochemistry - Plastic Embedding and Staining (Sharon Van de Velde, *B.S.*)
- Metabolic Bone Disease in Gastrointestinal Disorders (Dhanwade Rao, *M.D.*)
- Immunohistochemical Identification of Cells in Contact Dermatitis (Jules Ellis, *M.A.*)
- Undecalcified Bone Stains - Review and Comments of Certain Stains (Antonio Villanueva, *Ph.D.*)
- Discussion, Questions and Answers

Symposium Registration Form

DO NOT USE
THIS SPACE

Name: (last) (first) (initial)

Home Address: (street) (city) (state) (zip)

Employer: Address: (street) (city) (state) (zip)

Work Telephone No. () area code

Are you an NSH Member? Yes No

Is this your first attendance to an NSH Symposium/Convention? Yes No

On application below, please indicate first and second choice on limited workshops.

Check to ensure registration is for ONE, all day workshop, OR an AM and PM combination.

LATE REGISTRATION FEE OF \$10 required after September 7, 1982.

CANADIAN registrants please remit fees in U.S. currency.

Make check/money order payable to: National Society for Histotechnology

MAIL complete form & fees to: NSH, P.O. Box 36, Lanham, Maryland 20706

MONDAY			TUESDAY			WEDNESDAY-THURSDAY-FRIDAY		
	Member Fee	Non-Member Fee		Member Fee	Non-Member Fee		Member Fee	Non-Member Fee
All Day #1	\$40	\$45	All Day #14	\$40	\$45	Scientific Sessions	\$50	\$60
All Day #2	\$40	\$45	All Day #15	\$40	\$45	(Sessions & fee include Wed. thru Fri.)		
All Day #3	\$40	\$45	All Day #16	\$40	\$45	BANQUET \$15		
All Day #4	\$40	\$45	All Day #17	\$40	\$45	TOTAL REGISTRATION		
AM #5	\$20	\$25	AM #18	\$20	\$25	Monday Workshops	\$	
AM #6	\$20	\$25	AM #19	\$20	\$25	Tuesday Workshops		
AM #7	\$20	\$25	AM #20	\$20	\$25	Scientific Sessions		
AM #8	\$20	\$25	AM #21	\$20	\$25	Banquet		
AM #9	\$20	\$25	PM #22	\$20	\$25	Late Fee \$10 (required after Sept. 7, 1982)		
PM #10	\$20	\$25	PM #23	\$20	\$25	TOTAL FEES:		
PM #11	\$20	\$25	PM #24	\$20	\$25			
PM #12	\$20	\$25	PM #25	\$20	\$25			
PM #13	\$20	\$25						
Totals:								

SCIENTIFIC SESSIONS

WEDNESDAY, SEPTEMBER 15, 1982

HISTOTECHNOLOGY: THEN AND NOW
HANDLING HAZARDOUS MATERIAL IN THE LABORATORY
ESTROGEN RECEPTOR ASSAY IN BREAST CANCER
HT/HTL PRACTICAL EXAMINATION CRITERIA
HISTOCHEMICAL COMPARISON OF SENILE PLAQUES IN MAN
AND DOG

Israel Diamond, M.D.
Guy Fragala, Ph.D.
Neptune Mizrahi
William Kingsley, M.D.

Robert Pinzok, M.S.

PROFESSOR C.F.A. CULLING LECTURE

IMMUNOHISTOLOGY - A NEW FRONTIER IN SURGICAL
PATHOLOGY
CHRONIC AND LATENT EFFECTS OF SODIUM ARSENITE IN THE
AMERICAN OYSTER (*Crassostrea Virginica*)
WHAT EVERY HISTOTECHNOLOGIST SHOULD KNOW ABOUT
MALIGNANT MELANOMA
SUDDEN DEATH IN HEART ATTACK
NSH MEMBERSHIP MEETING (Opening Meeting) 4:45 P.M.

Clive Taylor, M.D.

Carolyn Barszcz, B.S., HT (ASCP)

John Kaiser, M.D.
Norman McLetchie, M.D., FRCPM, FCAP

THURSDAY, SEPTEMBER 16, 1982

CAP WORKLOAD RECORDING SYSTEM
THE USE OF HISTOPATHOLOGICAL TECHNIQUES IN THE
DIAGNOSIS OF FISH DISEASES
PROBLEMS IN DIAGNOSTIC HISTOPATHOLOGY - CASE STUDIES
ROLE OF HISTOLOGY LABORATORY IN FORENSIC MEDICINE
SLOW VIRAL INFECTIONS IN THE CENTRAL NERVOUS SYSTEM
MARINE PATHOLOGY
HISTOLOGY LABORATORY HAZARDS: FIRE AND FORMALDE-
HYDE
MORTALITY STUDY OF HISTOLOGIC TECHNICIANS CERTIFIED
BETWEEN 1948-1970 - FINAL RESULTS

Suzanne Carbone, HT (ASCP) & Thomas Webb

Arthur Hauck, M.S., PHI (AFS)
Earl Kasdon, M.D.
Luke Tedeschi, M.D.
William Pendlebury, M.D.
Paul Yevich, B.S.

Nicholas Hardin, M.D.

Dan Grauman, M.A.

FRIDAY, SEPTEMBER 17, 1982

SKIN IS IN - SPECIAL DERMATOPATHOLOGY TECHNIQUES AND
WHY EVERY HISTOTECH SHOULD KNOW THEM
METHODOLOGY AND INTERPRETATION OF FINE NEEDLE
ASPIRATION
UPDATE ON THE EFFECTS OF IN UTERO EXPOSURE TO
DIETHYLSTILBESTROL IN FEMALE AND MALE OFFSPRING
FROM SPILSBURY TO QUINCY: SCIENTIFIC DETECTION OF
CRIME - THE DIAGNOSIS OF MURDER
SPORTS MEDICINE
LECTURE from Paper Written During Workshop on Tuesday

Florence Nygaard

Rosa Enriquez, M.D.

Donald Antonioli, M.D.

William Sturner, M.D.
Arthur Pappas, M.D.

less than its TLV-TWA. For example, xylene's TLV-TWA is 100 ppm, but its ROTC is 0.47 ppm — about 200 times less than its TLV-TWA. Thus, while the concentration of xylene in laboratory air is usually less than 100 ppm and is ordinarily within "safe" levels, its presence is readily detected by smell at far less concentrations. Smelling even very low concentrations of xylene can cause headaches. In such instances, the concentration of xylene is unacceptable, even though within "safe" limits.

The vapor concentrations of formaldehyde and HCl, as well as other toxic chemicals, can be maintained below their TLVs-TWA by taking the following measures: (1) minimize their exposure to air by keeping all containers tightly capped or securely covered when the solutions are not in use; (2) use the solutions away from evaporation-promoting circumstances such as draughts from air conditioners, people traffic, proximity to doorways and sources of heat; and (3) use these chemicals only within the work space of a conventional, properly functioning fume hood.⁴

Self-contained partial enclosures with charcoal sorbents (e.g., Labconco's Fume Absorber,⁵ Lerner Laboratories' Fume-Gard,⁶ and Lipshaw's Fume-X-Peller⁷) cannot reduce the concentration of vapors as effectively as conventional fume hoods and are not recommended. Parenthetically, these units are misnamed, for it is not fumes they attempt to control, but vapors. Vapors are the gaseous forms of substances which are normally in the solid or liquid state (at room temperature and pressure). On the other hand, fumes are solid particles generated by condensation from the gaseous state, generally after volatilization from molten metals.⁸ Fumes are not ordinarily produced in histology laboratories.

MDA Scientific's Captair 2007C is also a self-contained partial enclosure designed to filter toxic vapors. It uses a catalytic-reaction molecular filter as a sorbent and is more effective than charcoal. Its filtration efficiency has been demonstrated by an independent testing laboratory to be nearly 100 percent.⁹ The Captair unit appears to be a worthwhile alternative to the expense of installation and operation of a conventional fume hood.

In summary, formaldehyde and hydrochloric acid are highly unlikely to yield bis-chloro-methyl ether under ordinary conditions of use in a histology laboratory. B-5 fixative is an improbable source of BCME. The concentration of a chemical's vapors in air should be kept so low that its presence escapes detection by smell. Self-contained recirculating partial enclosures that use charcoal sorbents to remove toxic vapors from the air are not recommended. Toxic chemicals such as formaldehyde and hydrochloric acid should be used under a conventional fume hood to ensure personnel safety.

References:

1. Frankel, L.S., McCallum, K.S., and Collier, L.: Formation of Bis (chloromethyl) Ether from Formaldehyde and Hydrogen Chloride. *Env. Science Technol.* 8:356-359, 1974.
2. American Conference of Governmental Industrial Hygienists: Threshold Limit Values for Chemical Substances in Workroom Air Adopted by ACGIH for 1975. ACGIH, Cincinnati, 1978.
3. Lillie, R.D., and Fullmer, H.M.: *Histopathologic Technique and Practical Histochemistry*, 4th Edition. McGraw-Hill, New York, 1978.
4. Scientific Apparatus Makers Association: Standard for Laboratory Fume Hoods, SAMA, Washington, D.C., 1975; 1980 revision available end of 1980. (Free copy available by writing SAMA, 1101 16th Street, N.W., Washington, D.C., 20036, or by calling (202) 223-1360.)
5. Labconco Corporation, 8811 Prospect, Kansas City, Mo. 64132.
6. Lerner Laboratories, 17 James Street, New Haven, Conn. 06513.
7. Lipshaw Manufacturing Corporation, 7446 Central Avenue, Detroit, Mich. 48210.
8. American National Standards Institute: Fundamentals Governing the Design and Operation of Local Exhaust Systems - ANSI Z9.2-1979. ANSI, New York, 1979.
9. MDA Scientific, Inc., 800 Busse Highway, Park Ridge, Ill. 60068.

Decalcifying Fluid Explosion

Leigh Winsor

James Cook University
School of Biological Sciences
Queensland 4811, Australia

The following is a recent report by Pirie and Soretire¹ regarding an explosion of a winchester of nitric acid-formalin decalcifying fluid:

"... an amber winchester half full of a 10% nitric acid-formalin decalcifying fluid exploded spontaneously with such force as to cover a room approximately 22 feet by 26 feet (6.6m x 7.8m) with fragments of broken glass. The bottle, which had been used earlier that morning, was standing in the place it had occupied for several weeks. Neither the room nor the weather was hot and the bottle was not in direct sunlight ... nobody was in the room at the time of the incident ... It is suggested that a build-up of gaseous methyl nitrite, formed by the reaction of traces of nitrous acid with the methanol usually added to prevent formaldehyde polymerization, is a possible cause."

Comment:

This incident bears certain similarities to explosions of sealed containers of chromic acid and aqua regia reported in the literature.² Possible explanations for these explosions include oxidation of contaminants and autodecomposition catalyzed by impurities present in the reagents involved. Similar mechanisms may, in general, apply to the nitric acid-formalin explosion.

Numerous decalcifying fluid mixtures and techniques have been thoroughly investigated, evaluated and rationalized in two excellent studies by Clayden³ and Brain.⁴ They have clearly demonstrated that 5-10% simple aqueous solutions of nitric acid or formic acid are the best choice of decalcifying fluids for routine diagnostic purposes. More complex and exotic mixtures were found to be generally inferior to these two simple fluids, providing that the latter were used correctly. It has also been found that substances, including formalin, added to nitric acid or formic acid decalcifying fluids to prevent maceration, serve no purpose when employed on well fixed tissues.⁴ Simultaneous fixation and decalcification of tissues saves the time required for separate fixation. However, this time-saving is more than offset by the inferior quality of the resulting tissue sections.⁴

Although more limited, my experiences with decalcifying fluids lead me to concur with the foregoing findings, and I would urge that laboratories which at present employ nitric acid or formic acid-formalin mixtures, replace them with the safe, single decalcifying fluids recommended. Should personal preferences and prejudices see the retention of the nitric acid-formalin type fluids, then safe use can be simply attained by preparing the mixture immediately before use, or by storing the decalcifying fluid in a vented container.

The report highlights the safety problem inherent in many of the classical brews which abound in Histotechnology. A trend towards the use of scientifically formulated, rationalized, simple and effective reagents would be welcome.

References:

1. Pirie, C.J., and Soretire, E.A.: Decalcifying Fluid Explosion. *Chemistry in Britain*, 15:11-12, 1979.
2. Winsor, L.: Potential Health Hazards and Safety in the Histology Laboratory: A Review. *Austral. J. of Med. Technol.* 9:149-166, 1978.
3. Clayden, E.C.: A Discussion on the Preparation of Bone Sections by the Paraffin Wax Method with Special Reference to the Control of Decalcification. *J. of Med. Lab. Tech.*, 10:103-123, 1952.
4. Brain, E.B.: *The Preparation of Decalcified Sections*. Charles C. Thomas, Springfield, Ill., 1966.
5. Wallington, E.A.: *Histological Methods for Bone*. Butterworths, 1972.



LAB-TEK DIVISION
 Miles Laboratories, Inc.
 30W475 North Aurora Rd.
 Naperville, IL 60566
 Address Correction Requested

NAME OR ADDRESS CHANGE? NEW SUBSCRIPTION?
 Please PRINT correct information below and return with old label (except new subscribers) to Lab-Tek. Mailing list updated quarterly. Allow 90 days for change to be made.

Name _____

New Address _____

City, State _____

Zip _____

☐ Please take my name off the Histo-Logic mailing list.
 (Attach old label.)

BULK RATE
 U.S. POSTAGE
 PAID
 PERMIT NO. 4954
 Chicago, Illinois

FREE INTRODUCTORY OFFER!

Buy the new Tissue-Tek® III Cassette Organizer now and get free Uni-Cassettes™




Supply is limited. Order now. Offer good only in U.S.A.

For a limited time only, we'll give you a dispenser box of 250 Tissue-Tek® III Uni-Cassettes (4170) with each purchase of a Tissue-Tek III Cassette Organizer (4136).

The new Tissue-Tek III Cassette Organizer makes it easy to label and store processing/embedding cassettes. Cassettes are held upright for convenient identification. The Cassette Organizer's front compartment provides storage for Tissue-Tek II and III Cassettes, process covers, marking pencils and small accessories. Rear compartment accommodates three Cassette Dispensers (4137, included). These dispensers hold either Tissue-Tek II or Tissue-Tek III Cassettes. Or use rear compartment to hold three Mega-Cassette™ Dispensers, available separately (4139). Cassette Organizer has a brushed aluminum lid and metal body; 11"x9"x3 1/4".

Cassette Organizer and Uni-Cassettes . . . the perfect partners for permanent, complete specimen identification from grossing to filing. Special offer in effect only as long as supply lasts. Offer good only in U.S.A.

So contact your distributor today and order the new Tissue-Tek III Cassette Organizer. Available through American Scientific Products, Curtin Matheson Scientific and VWR Scientific.



The histology systems source

To receive your own personal copy of HISTO-LOGIC, or to have an associate added to the mailing list, submit home address to: Lab-Tek Division, Miles Laboratories, Inc., 30W475 North Aurora Rd., Naperville, Illinois 60566

Printed in U.S.A.

The editor wishes to solicit information, questions, and articles relating to histotechnology. Submit these to: Lee G. Luna, Editor, Histo-Logic, P.O. Box 36, Lanham, Maryland 20706. Articles, photographs, etc., will not be returned unless requested in writing when they are submitted.