

Lee G. Luna, D. Lit., H.T. (ASCP), Editor

A Technical Bulletin for Histotechnology

Vol. V. No. 4 - October, 1975

Four-Hour Processing Schedule for Mouse Tissue

Wilson A. Werely, B.S., HT (ASCP) **Pathology Services** National Center for Toxicological Research Jefferson, Arkansas 72079

Our records show that in 1974 we necropsied, processed, embedded, sectioned and stained H&E slides on 98,847 adult and fetal mice. To handle this large volume of tissue it became necessary to develop a modified 4-hour processing schedule for the 10 Ultra tissue processors, used in the Histopathology Laboratory Branch.

The delay timer was modified to provide an overnight run of 450 Tissue-Tek® II Cassettes* (50 adult mice) and made ready for embedding in the morning when the histotechnologist reports for work. In addition, we scheduled a day run with 450 cassettes which are ready for embedding by 1:00 p.m. that same day. One to 7 organs are placed in each cassette and I to 10 cassettes per animal depending upon the experiment being conducted.

On the day schedule the various solutions are rotated and/ or changed. The following solution changes are made just prior to starting the day processing schedule. (A) Replace the 70% and 80% alcohols (beakers 1 and 2) with fresh alcohols. (B) Discard the 95% and absolute alcohols (beakers 3 and 5), (C) Move 95% alcohol (beaker 4) to beaker station 3, Replace beaker station 4 with fresh 95% alcohol. The same rotational procedure is used with the absolute alcohol. Move beaker from stations 6 and 7 to beaker stations 5 and 6. Replace beaker 7 with fresh absolute alcohol. The same procedure can be followed for the xylenes. (D) The impregnating paraffins are replaced completely for both the day and overnight schedules.

The overnight run is loaded on the Ultra in the afternoon and the delayed timer used. For this run all alcohols and infiltrating media are replaced. To properly infiltrate mouse tissue, infiltration media must be kept fresh.

The processing temperature used is 35" C, paraffin baths are 57" C to 60" C. Animal tissue should not be exposed to temperatures over 60" C.

	Four-Hour Processing Sch	hedule
Beaker	Solutions	Time
1	70% ethyl alcohol	25 minutes
2	80% ethyl alcohol	10 minutes
3	95% ethyl alcohol	10 minutes
4 5	95% ethyl alcohol	10 minutes
5	100% ethyl alcohol	15 minutes
6 7 8 9	100% ethyl alcohol	20 minutes
7	100% ethyl alcohol	20 minutes
8	Xylene	10 minutes
9	Xylene	15 minutes
10	Xylene	15 minutes
11	Paraffin	20 minutes
12	Paraffin and Paraplast (2 parts) (1 part)	50 minutes

*Lab-Tek Products, Naperville, Illinois

Mention of products in this paper does not constitute an endorsement by NCTR, FDA, or EPA.

COPYRIGHT (11979, LAM. TEX PRODUCTS, DIVSION MILES LABORATORIES, INC.

Automated Method for Deceration and Rehydration of Microscopic Slides

Jean Williams and Rena S. Futch Florida Department of Natural Resources Marine Research Laboratory St. Petersburg, Florida 33701

Although there are automatic staining machines on the market today, it may be impractical to employ them in small clinical or research laboratories because of limited work loads. Fortunately, many of these laboratories have automatic tissue processors which can be used to perform a dual function,

Our laboratory uses a double decker Autotechnicon in the following manner: The lower level is set up for tissue processing while the upper level is used for decerating and rehydrating slides according to the following schedule.

- 1. Xylene, 5 minutes.
- 2. Xylene, 5 minutes.
- 3. Xylene, 5 minutes.
- 4. Absolute alcohol, 2 minutes.
- 5. 95% ethanol, 2 minutes. (See remarks)
- 6. 70% ethanol, 2 minutes. (See remarks)
- 70% ethanol, 2 minutes,
- 8. Distilled water, 2 minutes.

Slides are transferred to coplin jars or staining dishes filled with distilled water or 70% ethanol, depending on which stain is to be used.

Remarks For Bouin's fixed material, insert a beaker of saturated sodium thiosulfate in 70% alcohol between steps 5 and 6. This solution removes the picric acid from the slides and insures better staining reactions. Exposure time to this solution should be 5 minutes.

If a mordant is required, as is the case with some special stains, that solution is inserted between steps 6 and 7. The exposure time is dictated by the stain used.

The different timing schedules are determined by notching the timing discs according to procedure used. It only takes a few minutes to reposition the beakers and set the timing disc. Slides which are being stained in a variety of stains are identified by using the stain abbreviations or number coded. A single level tissue processor can be used for deceration and rehydration if one maintains an extra set of beakers filled with the proper solutions, These beakers replace the processing beakers and the properly set timing disc is used.

Decalcifying Solution for Hard or Soft Tissue'

Livia M. Molnar Department of Orthodontics University of Washington Seattle, Washington 98195

Decalcification is thorough and rapid with this solution (1-7 days depending on the size of the hard or soft tissues) and there is no need to change the solution every day. Solution can

be kept in a dark bottle on a shelf for six months in good standing.

		olution

Distilled water	÷	A.	4	-	.,		2	4	÷,	. 1	Im 0.0101
Formic acid .											90.0 ml
HC1											80.0 ml
Sodium citrate											10.0 gm

Decalcifying Procedure

 Fix tissue specimens in any type fixative. (Formalin preferred)

- 2. Transfer tissue to decalcifying solution.
- 3. Let stay for 1-7 days (depending on size of specimen).
- Check with x-ray every 3rd day (again, depending on size of specimen).

Preparation of this note was supported in part by grant RR00166 from the National Institutes of Health.

¹We wish to acknowledge Dr. K. Bankuthy (deceased) who developed this solution and shared it with our laboratory

Additional Hints

Brenda Cuevas and Jose Torres Gorgas Hospital Laboratory Balboa Heights, Canal Zone

We have read with great interest the hints which appeared in the AFIP Letter,* Vol. 133, No. 9, December 1974, paragraph four. We would like to share our procedure with readers of HISTO-LOGIC. The following solution is used in this procedure:

Glycerin-Alcohol Solution

Alcohol,	70%	 90.0 ml
Glycerin		 10.0 ml

Softening Procedure

- Trim (rough cut) the paraffin block until the tissue is exposed.
- Place paraffin block in a petri dish containing glycerinalcohol solution and let stand for 1 hour.
- Remove paraffin block from solution and remount on the microtome.
- 4. Chill block surface with a piece of ice and section.

* The AFIP Letter referred to contained the following information: Dried tissue, thyroid colloid, hemorrhage, eye lens, bone, etc., can be sectioned intact if some water is applied to the tissue in the following manner: Rough cut block until all tissue is exposed. Dip a piece of cotton into the flotation bath and squeeze out excess water. Apply water-soaked cotton to the surface of the block for a few seconds. Chill paraffin block surface with a piece of ice and section.

ALSO

In HISTO-LOGIC, Vol. 4, No. 2, April 1974, Shirley Orlando requested ideas pertaining to the difficulty of obtaining sections from specimens containing large amounts of fat; lipomas for example. We have developed the following procedure and find it most useful in dealing with this type of problem.

Solutions

	10% Alcoholic Formalin	
For	maldehyde, concentrated (37-40%)	10.0 ml
Alc	ohol, 80%	90.0 ml
	Acetone-Chloroform	
Act	tone	100.0 ml
Chl	oroform	100.0 ml

Procedure

- 1. Tissue specimen should not exceed 2 mm in thickness.
- Fix tissue in 10% alcoholic formalin solution for a minimum of 6 hours and preferably 24 hours.
- Place tissue in acetone-chloroform solution for 2 hours in a 60th C oven.
- Remove tissue from this solution and process in the usual manner.

Modified H&E Stain for Demonstrating Reed-Sternberg Cells

Peter V. Emanuele Histopathology Laboratories Division Armed Forces Institute of Pathology Washington, D.C. 20306

This modified hematoxylin and cosin staining procedure has been used to good advantage in the demonstration of acidophilic nuclear material present in Reed-Sternberg cells.

Fixation Any well fixed tissue (10% buffered neutral formalin preferred).

Microtomy Cut sections at 3-4 microns

Solutions

Mayer's Hematoxylin

Hematoxylin crystals	1.0 gm
Distilled water	1000.0 ml
Sodium iodate	0.2 gm
Ammonium or potassium alum	50.0 gm
Citric acid	1.0 gm
Chloral hydrate	50.0 gm

Dissolve the hematoxylin in the distilled water, add the remaining chemicals in the order given, making sure each is completely dissolved before the next one is added. Use a magnetic stirrer. The final color of the stain is reddish-violet. Stain keeps well for months.

1% Phloxine B (stock)

Phloxine	Β.	 	 4.14	-	 + +	+ +	+	÷.,	 1.0 gm
									100.0 ml

Eosin-Phloxine B Solution (working)

1% Eosin (stock)			+	100.0 ml
1% Phloxine B (stock)	-	+ + +	•	10.0 ml
Alcohol, 95%			ł	780.0 ml
Glacial acetic acid,	-		-	5.0 ml

Make up working solution as needed. Working solution should be changed at least once a week.

Staining Procedure

- 1. Decerate and hydrate to water in the usual manner.
- If sections are Zenker-fixed, remove the mercuric chloride crystals with iodine and clear with sodium thiosulfate (hypo) in the usual manner.
- 3. Stain in Mayer's hematoxylin for 15 minutes,
- 4. Wash in running tap water for 20 minutes.
- Counterstain with eosin-phloxine B for 8 minutes. For even staining, dip slides several times before allowing them to sit in the eosin for the prescribed time.
- Dehydrate in 95% and absolute alcohols, 2 changes each. Pass slides through these solutions as rapidly as possible, making sure all water is removed.
- 7. Clear in xylene, 3 changes of 2 minutes each.
- 8. Mount coverglass in resinous media.

Results Nuclei – blue Reed-Sternberg cell acidophilic nuclear material – bright red Cytoplasm – various shades of pink

National Society for Histotechnology Symposium/Convention

The Second Annual Symposium/Convention of the National Society for Histotechnology will be conducted at the Sheraton and Holiday Inn. Silver Spring, Maryland, October 6-10, 1975.

Registration application is attached. We request filing your egistration by September 15th, as the scientific sessions are amited to 600 registrants. Late applications will be accepted provided quota has not been filled. Registration forms may be photocopied if more than one individual from the same activity wishes to attend. Mail registration and check to: Registrar, P.O. Box 36, Lanham, Maryland 20801.

Also attached are hotel reservation cards for the Sheraton and Holiday Inn. Please complete reservation card of choice and mail with registration to: Registrar, P.O. Box 36, Lanham, Maryland 20801, for recording purposes. As our workshop

room rentals are determined by the number of registrants staying in the hotel, this will enable the NSH to have better negotiating power for prices. Reservation card will then be submitted to respective hotel who will send confirmation directly to you.

SHERATON RESERVATION CARD MUST BE USED as their Reservatron machine will be closed for this meeting. RESERVATION WILL NOT BE ACCEPTED IF TOLL FREE RESERVATRON IS USED. You may call Sheraton directly for late reservations.

Listed below are the titles for workshops, scientific sessions, and seminars being conducted:

(4) Common problems and solutions. (5) A question and

Workshops

Topics for discussion during this workshop include, Communications in Management, Perception and the Behavioral

Sciences, Listening as a Factor in Communications, How Can

I Be a Better Manager, and a Workshop Summary, Learn to be

a better manager by discussing and employing the latest behav-

ioral sciences techniques with your peers. Each workshop

9:00 A.M. - 4:30 P.M. (Sheraton)

MANAGEMENT WORKSHOP Dr. James Young

pate in the management process.

WORKSHOP NO. 5

answer period.

TUESDAY, OCTOBER 7, 1975

9:00 A.M. - 4:00 P.M. (Sheraton)

CRYOSTAT TECHNOLOGY

Mr. Frank Avallone

This workshop is designed to include information on basic and advanced techniques in the field of cryotomy, Lectures will be presented on different aspects of cryotomy and include such topics as general and specialized techniques in the field of histochemistry and immunofluorescence, routine staining procedures, and a short history of cryotomy. There will be a problem session and a limited number of cryostats for class use, Registrants are encouraged to discuss their problems, whether they deal with knives, tissues, or some other aspect of cryostat technology, during the workshop meeting.

WORKSHOP NO. 6

TUESDAY, OCTOBER 7, 1975

TUESDAY, OCTOBER 7, 1975

9:00 A.M. - 12:00 NOON (Sheraton)

FLUORESCENT ANTIBODY TECHNIQUES Mr. C. F. A. Culling

(No abstract received by date of printing.)

WORKSHOP NO. 7

9:00 A.M. - 12:00 NOON (Sheraton)

MANAGEMENT IN HISTOLOGY LABORATORY Ms. Elizabeth Mayle

Histotechnology laboratories in certain areas can be treated like any other organized laboratory setting. In many areas there are peculiarities and differences that need to be defined, investigated and improved upon, so that we may have more professionalism and efficiency.

This workshop will define particular difficult areas for both the one-technologist laboratory and the multi-technologist laboratory, whether it be research oriented or patient care oriented. It will give specific suggestions and guidelines to the histotechnologist so that these difficulties can be alleviated.

WORKSHOP NO. 8 **TUESDAY, OCTOBER 7, 1975**

1:00 P.M. - 4:00 P.M. (Sheraton)

PROCEDURES FOR DETERMINATION OF DRUGS IN BLOOD AND URINE

Dr. Maurice Goldbaum (No abstract received by date of printing.)

WORKSHOP NO. 4

TUESDAY, OCTOBER 7, 1975

9:00 A.M. - 4:00 P.M. (Holiday Inn)

SPECIAL STAIN SEMINAR Mr. Jack Wenger

This seminar on stain technology will deal with: (1) A review of staining theory including reactions of silver and other methods. (2) Preferred stains for demonstrating specific entities. (3) Permissible variables and innovative stain technology.

9:00 A.M. - 4:00 P.M.

(Sheraton)

MICROTOMY

Mr. H. Keith Russell

The problems relative to poor versus good microtomy will

be emphasized. The various steps in cutting, which may be responsible for variable results, will be discussed, including the effects on tissue sectioning of the following: Fixation, processing, embedding, microtomy, and tissue orientation. Indepth knife sharpening methods will be included. Participants are encouraged to bring "problem" blocks and microtome blades, as you will be sharpening blades, cutting tissue, staining slides, and taking home the finished product.

member will be provided the opportunity to actively partici-

MONDAY, OCTOBER 6, 1975

9:00 A.M. - 4:00 P.M. (Holiday Inn)

INTRODUCTION TO TISSUE IDENTIFICATION Mrs. Edna Prophet and Mr. Charles Collard

The purpose of this workshop is to acquaint the participant with basic tissue structures. The value of this knowledge in microtomy and staining technology will be illustrated with the use of photomicrographs.

WORKSHOP NO. 3

VORKSHOP NO. 2

MONDAY, OCTOBER 6, 1975

WORKSHOP NO. 9 TUI

TUESDAY, OCTOBER 7, 1975

8:30 A.M. - 11:30 A.M. 1:30 P.M. - 4:30 P.M.

(20 Histotechnologists)

(20 Histotechnologists)

SPECIAL STAIN WORKSHOP (Sheraton) Mr. Erwin Haas

This workshop will cover the following special stains: Verhoeff's Elastic; Gordon Sweet Reticulum; Masson Trichrome; Benhold's Congo Red; McCallum-Goodpasture's; Ziehl Neelsen. The workshop will be given twice a day, from 8:30 to 11:30 A.M. for twenty Histotechnologists, and from 1:30 to 4:30 P.M. for twenty Histotechnologists. If necessary, the workshop will be conducted on Monday and Tuesday to accommodate 80 Histotechnologists. A question and answer period and slide presentation will follow the practical exercise for Histotechnologists wishing a more in-depth discussion of the special stains performed during the practical exercise.

WORKSHOP NO. 10 TUESDAY, OCTOBER 7, 1975

6:30 P.M. - 7:30 P.M. (Sheraton)

SEMINAR ON FORMATION OF STATE HISTOTECHNOLO-GY SOCIETIES

Mr. Donald Hammer

Seminar will present methods, hints and procedures helpful in organizing a State Histotechnology Society. Individuals interested in obtaining information on forming a Society in their state should avail themselves of this opportunity.

Scientific Sessions

(All Sessions Held At The Sheraton)

WEDNESDAY, OCTOBER 8, 1975:

AN ALTERNATE APPROACH TO STAINING ACID FAST ORGANISMS

Mr. Wayne Kampa

THE ROLE OF HISTOPATHOLOGY IN THE MEDICO-LEGAL INVESTIGATION Dr. Richard Froede

QUALITY CONTROL IN THE HISTOPATHOLOGY LABORATORY

Ms. Elaine Boyd

CATECHOLAMINE HISTOFLUORESCENCE Mr. Robert Welsh

TELLTALE TISSUE TRACES Ms. Barbara Campbell

IMMUNOFLUORESCENT VISUALIZATION OF GASTRO-INTESTINAL ENDOCRINE CELLS, USING PARAFFIN EM-BEDDED BIOPSIES Mr. Robert Escoffery

HISTO- AND CYTO-TECHNOLOGY IN VETERINARY MED-ICINE

Mr. Thomas Palmer

WEDNESDAY EVENING:

7:00 - 9:00 P.M. NATIONAL SOCIETY FOR HISTOTECH-NOLOGY BUSINESS MEETING, EVERYONE WELCOME!

THURSDAY, OCTOBER 9, 1975:

HISTOTECHNOLOGY IN AN INDUSTRIAL PATHOBIOL OGY OPERATION – A UNIQUE APPROACH Dr. Gary Johnson EXPERIENCES WITH AMYLOID STAINS AND RECOM-MENDATIONS FOR IMPROVING THEIR USE AND INTER-PRETATION Dr. Robert Mowry

THE RUSSELL/MOVAT STAIN, ITS VALUE TO THE

PATHOLOGIST Dr. Daniel Connor

THE AMA AND THE ACCREDITATION PROCESS Dr. John Beckley

MUSCLE ENZYME HISTOCHEMISTRY Ms. Alice Harris

CONTROL OF AND TROUBLESHOOTING FOR SPECIAL STAINING PROBLEMS Mr. John Koski

SILVER METHODS FOR RETICULIN AND CARBOHY-DRATES ARE SPECIFIC, AREN'T THEY? Mr. C. F. A. Culling

THURSDAY EVENING:

6:30 – 7:30 P.M. Complimentary cocktails sponsored by LAB-TEK PRODUCTS, Naperville, Illinois

7:30 - 9:30 P.M. National Symposium Banquet

Banquet Speaker: Ralph C. Kuhli, MPH, Director, Department of Allied Medical Professions & Services, American Medical Association, Chicago, Illinois.

Subject Title: "NATIONAL COOPERATION FOR ALLIED MEDICAL EDUCATION"

FRIDAY, OCTOBER 10, 1975:

SUCCESSFUL LOBBYING – HOW TO CONVINCE YOUR CONGRESSMAN, HOSPITAL, PATHOLOGIST, TO DO WHAT YOU WANT Mr. Keith Russell

ALDEHYDE FUCHSIN STAINING OF ELASTIN AND MU-CINS

Ms. Beatrice Macdonald

A PERIODIC ACID-METHENAMINE SILVER-CHROMO-TROPE 2R METHOD FOR STAINING BASEMENT MEM-BRANE

Mr. Wayne Kampa

PREPARATION OF OSSEOUS TISSUE EMPLOYING THE WET CELLOIDIN TECHNIQUE Mr. Gerald Armstrong

PROBLEMS OF PSEUDO-MINERALIZED BRAIN Mr. Celester Carter

HISTOCHEMISTRY Dr. Frank Johnson

HISTOLOGICAL TECHNIQUES USEFUL IN PARASITE IDENTIFICATION Mr. Ronald Neafie

AN EVALUATION OF THE PROBLEMS IN LYMPH NOE PREPARATION Mg. Elaine Boyd

CLOSING REMARKS Mr. Lee G. Luna

Day of above lectures subject to change pending finalization of program.

NAME_				I	O NOT USE THIS SPAC	E A
C. S.	Last	First	In	iitial		
HOME	ADDRESS					
	Street	City	State	Zip		
EMPLO	YER					
ADDRE	SS					
	Street	City	State	Zip	WORKSHO	
	CHIPPER PUNCTIONS NOT	INCOME TO ANTING		N	lo. 2\$20.00 (M	onday) onday) onday)
	CHECK FUNCTIONS YOU TFIC SESSIONS \$33			N	lo. 4 \$20.00 (T)	uesday)
	ET\$13.00 (Thurs			N	lo. 6\$20.00 (T)	uesday) uesday) ½ day uesday) ½ day
entrete	OR MONTH ORDER MINE			N	lo. 8 \$20.00 (T)	uesday) ½ day uesday AM
	OR MONEY ORDER MUST PAYABLE TO: NATIONAL				(T)	uesday PM
	EGISTRATION TO: REGI				lo. 10 No Fee (T	uesday Evening
		SHER. Colesville Road Silver Sp	ATON - SILVER		ana (101) 580 5200	
9 4	6121		OM RESERVATIO			
	NAME	UUEST KO	OM RESERVAIN	ON REQUES	•	
N	ADDRESS		a la constante			0.31
ENTI 1975	CITY	Section Sectio	CTATE.		ZIP_	
20	APRIL DATE					
R 6-	Reservations must be received OF ARRIVAL, unless guarant	a not later than two weeks p reed in writing. CHECKOUT	TIME IS 1 P.M.	nd ROOMS W	ILL BE HELD ONLY UNTIL	or P.M. on DATE
OCTOBER 6-1	ARRIVAL DATE			но	URA.M	P.M.
0410	DEPARTURE DATE		2.2.4.4	но	UR A.M	P.M.
SYI	SINGLE\$21.00	DOUBLE\$26.00	Suite by Special	Request Onl	y.	
	(All rates plus 7% Marylan	d Tax)				
	PLEASE NOTE THAT MO	OST SYMPOSIUM ACTIV	VITIES WILL BE	TAKING PL.	ACE AT THE SHERATO	N.,
	MAIL RESERVATION TO	D: REGISTRAR, P.O. B	OX 36, LANHAM	, MARYLAN	ID 20801	
	****************		**********			
			NN* - SILVER S			
GUM/CONVENTION OBER 6-10, 1975	8777	Georgia Avenue Silver S				
		RES	SERVATION REQ	UEST		
NO	NAME					
SYMPOSIUM/CONVENTION OCTOBER 6-10, 1975	ADDRESS					
10.1	стту	Section 2 -	STATE		ZIP_	
R 6-1	Reservations will be held unti	6 P.M., unless guaranteed o	r prepaid. Checkout (time is 12 Noo	n. Indicate type of reservatio	n below:
IPOSIUM/C OCTOBER	Guaranteed reservations will Prepaid reservations will	s will be held all night and g	uest will be billed for	full amount o	f room if he fails to show up.	
PCT 0	ARRIVAL DATE	and the second second		TIM	E AM	P.M.
o o	DEPARTURE DATE					
	SINGLE518.00					
	*Holiday Inn approximate	ly three blocks from She	raton, meeting her	adquarters		
		- Co.				

Russell's Fixative Solution

An Editorial

The substitution of zinc chloride for mercuric chloride in Zenker's fluid.

Distilled water	1000.0 ml
Zinc chloride	50,0 gm
Potassium dichromate	25.0 gm
Sodium sulfate	10.0 gm
Acetic acid, glacial	50.0 ml

This modified Zenker fixative utilizes zinc chloride in place of mercuric chloride. The use of mercuric chloride in fixative solutions should be discontinued due to the potential toxic effects of this chemical to the environment. The staining results produced by this fixative are similar to those produced by Zenker fixative containing mercuric chloride. Cadmium chloride may also be substituted for zinc chloride. Zinc and cadmium chloride are removed from sections by the usual dehydrating procedures; that is, no extra steps are required to remove crystals of the salts as are required after fixation with Zenker's solution containing mercuric chloride.

It is interesting to note that this procedure was developed in 1941 because of the high cost of mercuric chloride -- not the environmental problem posed by mercuric chloride.

Reference

Russell, W. O.: J. Tech. Meth. Bull., Int. Assoc. Med. Mus., 21:47-49, 1941.

-

Mounting and Storing Frozen Tissue Specimens

Marie I. Doman Department of Oral Biology University of Washington Scattle, Washington 98195

The following method is used in our laboratory for mounting and storing tissue specimens destined for cryotomy. The method has proved most useful and convenient.

Method

- Collected tissue specimens are oriented in O.C.T.* embedding media on a 1/8 inch slice of cork and plunged into liquid nitrogen.
- The frozen specimen can now be stored in the freezing compartment of the refrigerator in an air tight container or mounted on the cryostat metal object holder.
- If the tissue requires sectioning all the way through the block, there is no risk of damaging the knife edge on the metal object holder.
- 4. The specimen cut surface can be recovered with O.C.T., the cork removed from the object holder and the mounted specimen stored in a covered jar. If additional sections are required at a later date, the specimens will still be in their original orientation when the cork is remounted on the object holder.

This is done by first placing a drop of water on the object holder, then placing the specimen on the water for attachment to the object holder. The specimen can now be sectioned or stored while it remains attached to the cryostat object holder.

We find the following advantage to this method: It is a fast, handy method for orienting and storing large numbers of tissue specimens destined for cryostat technology.

We wish to thank Dr. James A. Claggett for his encouragement in writing this article.

*O.C.T. Compound, Lab-Tek Products, Naperville, Illinois

A Modified Argentaffin Silver Reaction

Karole Cecich Holy Family Hospital Des Plaines, Illinois 60016

The following silver nitrate procedure for argentaffin granules and melanin pigment is a modification which provides the following advantages:

- The method can be performed faster than other argentaffin type procedures.
- 2. It is simpler to perform.
- 3. The method gives excellent results.

This method has been used for the demonstration of melanin pigment, argentaffin and chromaffin granules.

Fixation 10% buffered neutral formalin

Microtomy Cut sections at 10 micra.

Solutions

Ammoniacal Silver Solution

To 10.0 ml of a 10% aqueous solution of silver nitrate add 2.0 ml of a 10% aqueous solution of potassium hydroxide. Add ammonium hydroxide (drop by drop) until all silver nitrate granules are dissolved. To this solution add 5.0 ml of a 10% silver nitrate solution. Pour above solution into 55.0 ml of distilled water. Filter solution before use.

10% Silver Nitrate

Silver ni	trate		 	 ****	10.0 gm
Distilled	water.	2.	 	 	100.0 ml

10% Potassium Hydroxide

Potassium hydroxide.	+ + +	 + + + + +	 10.0 gm
Distilled water	Ba	 1919-0-	100.0 ml

10% Formalin

5% Sodium Thiosulfate (Hypo)

Sodium	thiosul	fate	 	++++	*****	5.0 gm
Distilled	water	2.2.2				100.0 ml

Staining Procedure

- 1. Decerate and hydrate sections to distilled water.
- 2. Wash sections in running tap water for 5 minutes.
- 3. Rinse slides in 3 changes of distilled water.
- Impregnate slides in ammoniacal silver solution for 1 hour in a 58° C oven. Control silver deposition microscopically. Remove staining dish from oven when granules appear light brown.
- 5. Rinse slides in 3 changes of distilled water.
- Reduce slides in 10% formalin solution for one quick dip. This solution intensifies the reaction and turns the granules black.
- 7. Rinse slides in 3 changes of distilled water.
- 8. Place slides in 5% sodium thiosulfate for 5 minutes.
- 9. Rinse slides in 3 changes of distilled water.
- 10. Counterstain slides with cosin working solution of choice.
- Dehydrate in 3 changes of 95% alcohol, 2 changes of 100% alcohol, 2 minutes each.
- 12. Clear in xylene, 3 changes, 2 minutes each.
- 13. Mount coverslip with resinous media.

Results Argentaffin granules – black Background – pink

Remarks All glassware must be clean and free from anything which may contaminate the silver nitrate solution. Tissue has a tendency to fall off slides in ammoniacal silver solutions. A preventive measure to this problem is to use gelatin in the flotation bath and egg albumen adhesive on the glass slide. Chicago, Illinois

PERMIT NO. 4954 PAID

BULK RATE

A Technical Bulletin for Histotechnology Vol. V, No. 4



Address Correction Requested

Lab-Tek Products International Control of the Products

Announcing TISSUE-TEK® II TISSUE PROCESSOR



A hardworking processor to carry the whole load or an economical backup in an overloaded laboratory. Table-top size fits under exhaust hoods. Smooth-surfaced deck is easy to clean. All plugs are explosion-proof. Timing dial regulates entire cycle. Delay timer offers weekend economication up to 58 hours.

weekend processing up to 56 hours. Stainless steel perforated tissue baskets are in constant motion for thorough tissue penetration. Each holds up to

50 cassettes. Use in tandem for heavy loads. Insulated stainless steel paraffin pots have separate

thermostats.

Call your Lab-Tek dealer for a demonstration, or ask us directly.

Lab-Tek Products...setting the standards by which performance is judged.





HISTO-LOGICTM - Statement of Policy

There has been some misunderstanding by many people concerning who is eligible to receive HISTO-LOGIC. HISTO-LOGIC is available free of charge to any histotechnologist who wishes to receive it. You do not have to be a member of any organization or society, or falfill any prerequisites. The only requirement is that you make a request to be put on HISTO-LOGIC mailing list and farnish your name and forme address. Once you have been added to the mailing list, you will receive HISTO-LOGIC four times a year; January, April, July and October.

In addition, you will receive "NSH In Action." NSH In Action is the official newsletter for the National Society for Histotechnology. This newsletter contains news pertaining to all activities concerning the endeavors of the Society in the field of Histotechnology.

We ask that you encourage all histotechnologists to be placed on our mailing list. We also ask that you send all inquiries pertaining to HISTO-LOGIC to the Editor, P. O. Box 36, Lanham, Maryland 20801.

The editor feels that as Histotechnologists we are fortunate that Lab-Tek Products publishes HISTO-LOGIC and NSH in Action at no cost to the individual histotechnologist. I am sure most of you recognize this is a coaffy endeavor for Lab-Tek in both money and time.

The editor wishes to solicit information, questions, and articles relating to histotechnology. Submit these to: Lee G. Luna, Editor, Histo-Logic, P.O. Box 36, Lamham, Maryland 20801. Articles, photographs, etc., will not be returned unless requested in writing when they are submitted. 74

To receive your own personal copy of HISTO-LOGIC, or to have an associate added to the mailing list, submit home address to: Lab-Tek Products, Division Miles Laboratories, Inc., 30W475 N. Aurora Rd., Naperville, Illinois 60540.