

# HISTO-LOGIC<sup>®</sup>

No reader should utilize materials and/or undertake procedures discussed in HISTO-LOGIC articles unless the reader, by reason of education, training and experience, has a complete understanding of the chemical and physical properties of all materials to be utilized and the function of each material by which any procedure is accomplished.

Editor, Lee G. Luna, D. Litt., H.T. (ASCP)

Technical Bulletin for Histotechnology  
Published: January, April, July, October

Vol. XI, No. 3 - July, 1981

## Carbol-Xylene (A Useful Tool)

### An Editorial

Carbol\*-xylene was a very useful "tool" in the histopathology laboratory some years ago, primarily due to its softening effect on tissue sections embedded in celloidin. Sections obtained from celloidin-embedded specimens had a tendency to harden when exposed to xylene, resulting in section curling. This often produced wrinkles and, of course, distortion of tissue section. Carbol-xylene was used extensively during the celloidin era to soften these sections and allow them to uncurl and lie flat against the slide prior to the coverslip mounting. The strong emphasis on the use of paraffin over the last fifteen years has relegated carbol-xylene to the disuse status. However, this solution can still be very useful on occasion when one is experiencing detachment of sections from the slide, causing section "foldover." This is a very useful solution for tissues such as bone, chitin, or sections consisting of fibrous components, all of which have a tendency to dislodge from the glass slide.

### Solution:

#### Carbol-Xylene

Phenol (carbolic acid) ..... 25.0 gm  
Xylene ..... 75.0 ml

The procedure is simple to use by employing the following directions:

1. A camel's hair brush is dipped into the solution of carbol-xylene.
2. The brush is then applied to the folded portion of the specimen. (The application of carbol-xylene softens the section and allows easy unfolding by a gentle touch with the brush.)
3. Remove excess carbol-xylene from section. (If carbol-xylene is left on the section, tissue will fade.)
4. Mount coverslip with a resinous media.

\*Phenol or carbolic acid

## Liver Trichrome Stain

William Dotson

The North Carolina Memorial Hospital  
Chapel Hill, North Carolina 27514

The "liver trichrome" procedure is one of a battery of stains we do on all our liver biopsies. The pathologists seem to like it better than the Masson Trichrome, in this case, for its sharper color contrast. It is being presented here in the hope others will find it beneficial.

### Fixation:

10% buffered formalin or Helly's solution

### Microtomy:

Cut paraffin sections at 2 microns or desired thickness.

### Solutions:

#### Alcoholic Picric Acid (Stock)

Picric acid ..... 7.0 gm  
Ethanol, 95% ..... 100.0 ml

#### Alcoholic Picric Acid (Working)

Alcoholic Picric acid (stock) ..... 2 parts  
Ethanol, 95% ..... 1 part

#### 5% Aniline Blue

Aniline blue ..... 5.0 gm  
Acetic acid ..... 2.5 ml  
Distilled water ..... 100.0 ml

#### 1% Ponceau de Xylidine (Ponceau 2R)

Ponceau de xylidine (Harleco) ..... 1.0 gm  
Acetic acid ..... 1.0 ml  
Distilled water ..... 100.0 ml

#### 1% Acetic Acid

Acetic acid ..... 1.0 ml  
Distilled water ..... 100.0 ml

#### 1% Phosphomolybdic Acid

Phosphomolybdic acid ..... 1.0 gm  
Distilled water ..... 100.0 ml

### Staining Procedure:

1. Deparaffinize slides and hydrate to distilled water.
2. Stain slides in Harris' hematoxylin in 85°C for 5 minutes.
3. Dip slide once in 95% ethanol.
4. Place slides in working alcoholic picric acid for 10 minutes.
5. Wash slides in running tap water for 10 minutes.
6. Stain slides in ponceau de xylidine for 5 minutes.
7. Treat slides with 1% phosphomolybdic acid for 5 minutes (use solution only once).
8. Stain slides in 5% aniline blue for 2 minutes.
9. Rinse slides quickly in tap water to remove excess aniline blue.
10. Treat slides with 1% phosphomolybdic acid for 5 minutes.
11. Place slides in 1% acetic acid for 5 minutes.
12. Dehydrate, clear and mount coverslip with resinous mounting media.

### Results:

Connective tissue and collagen — blue  
Muscle and cytoplasm — red  
Nuclei — purple

## Making a Desk-Top Out of a Drawer

Jacob Lundy & Robert Moore

U.S. Army Biomedical Laboratory  
Aberdeen Proving Ground, Maryland 21010

It is frequently advantageous to be able to provide an extra working or writing surface at a counter or desk. Some desks have such a facility at one or both ends in the form of a pull-out top.

It is not difficult to make such a utility out of any drawer



under a counter top, especially in the modular-type cabinet systems now popular, because in these cabinets the drawer-wells are considerably larger than the drawers.




FIGURE 1.

A flat, rigid board, about 1 cm thick and as wide as the drawer, is obtained to rest upon the right and left sides of the drawer (Figure 1). A board at least three-quarters as deep (front-to-back) as the drawer provides an adequate working surface for many purposes, including use as a stand for a dual-viewing microscope. On the underside of the board, two small rubber feet are fastened (by screws or tacks) in a vertical line near the right and left sides, as guides (Figure 2). The rubber guides are positioned at a distance from the edges of the board to hardly touch the sides of the drawer with the board in place, thus assuring a free back-and-forth movement without sideways. The distance the board-top can be pushed back, and the area of inner drawer-space exposed, depends on the location of the guides at the rear of the board; that is, the top can be pushed back until the rear guides strike the back side of the drawer.

To determine the proper location of the 2 gear guides, the board is laid on the right and left sides of the drawer with its back edge coinciding with the edge of the back side of the drawer. If, with the board in this position, access to the drawer is adequate, the rear guides are attached against the rear edge on the underside of the board. If, however, a greater opening is desired, the rear guides are secured on the board at a distance from the rear edge equal to the additional millimeters of opening desired.




FIGURE 2.

A handy pull for manipulating the drawer-top can be made from a length of tape. Contact paper may provide a smoother surface and a more attractive appearance. In a cabinet where the drawer-well is too small to allow the entrance of a board atop the drawer sides, the board must be removed for reaching and looking into, and for closing, the drawer.

Note: The authors are thankful to Mr. Ray Clawson for the photography.



## Muscle Quenching with Liquid Nitrogen and Talcum Powder

Leonard Noble and Venkata Challa, M.D.  
North Carolina Baptist Hospitals  
Winston-Salem, North Carolina 27103

To achieve enzyme localization and to eliminate ice-crystal artifacts in muscle biopsies, the use of liquid nitrogen as a coolant is prescribed. Isopentane and various fluorocarbons can be cooled down by liquid nitrogen to a temperature which will adequately quench muscle biopsies. However, the use of these agents in freezing techniques has now become outdated due mainly to storage and health hazards.

We feel that the method of choice is to use liquid nitrogen alone. However, to accomplish this, something must be done to prevent the gaseous layer which forms around tissue when thrust into liquid nitrogen. This layer of gas acts as an insulator and does not allow for artifact-free freezing. For this reason, a coating of talcum powder is put on the tissue prior to freezing to eliminate the gaseous layer. As the powder falls away from the specimen, the gas layer is eliminated. The biopsy itself is anchored in place with a talcum powder-OCT paste which we have devised to take the place of the more often used gum tragacanth.

Moline and Glenner<sup>1</sup> reported the talcum powder coating principle for rapid quenching of tissue specimens in 1964. Their experiments clearly showed that adequate quenching of tissue is possible using liquid nitrogen alone and not in combination with hydrocarbons or fluorocarbons.

When adapted to muscle biopsies routinely, this technique is safe and relatively inexpensive. Freezing can actually take place in an area adjacent to the biopsy room, thus avoiding any contact with saline. After quenching, the specimen remains in the liquid nitrogen until it is transported to the laboratory.

### Materials:

Talcum powder  
Liquid nitrogen  
OCT embedding compound\*  
Cryostat chucks  
Cork discs (2.4 cm in circumference x 0.4 cm in thickness)†  
Dewar flask - 1 liter capacity

### Solutions:

#### Talcum powder-OCT paste

To make, mix small amount of OCT embedding compound with talcum powder until a paste-like consistency is obtained. Discard any leftover paste after freezing specimen.

### Method:

For best results, the biopsy should be 6 to 7 mm long and no more than 5 mm in cross section. If the dimensions are any longer, subsequent trimming may be necessary.

1. Cover one side of a cork disc with a generous amount of the talcum powder-OCT paste. Position the biopsy in the paste in such a way as to obtain a cross section of the muscle fibers when sectioning.
2. Coat the entire surface of the specimen as well as the paste and the cork disc with the talcum powder.



# National Society for Histotechnology Symposium/Convention

## November 16-20, 1981

## Salt Lake City, Utah

The Seventh Annual Symposium/Convention of the National Society for Histotechnology will be conducted at the Little America Hotel, Salt Lake City, Utah. The enclosed program is complete with hotel reservation card and registration form. The convention will utilize 425 sleeping rooms in the Little America, with overflow accommodations in the Tri-Arc Hotel, 801/521-7373 (which will also accommodate several workshops), and the Hilton Inn, 801/532-3344. Both hotels are within a half block of the Little America and are in walking distance. All room reservations will be processed through the Little America. When this hotel is filled, reservations will be forwarded to other hotels for accommodations. PLEASE MAKE YOUR RESERVATIONS EARLY SINCE ALL ROOMS BLOCKED FOR NSH WILL BE RELEASED ONE MONTH PRIOR TO MEETING DATE.

Mail Hotel Reservation Directly to: Little America, 500 South Main, Salt Lake City, Utah 84101; 801/363-6781.

Symposium registration application may be photocopied if more than one individual from the same activity desires to attend. To avoid delays and unnecessary complications, registrations AWAITING FUND APPROVAL will be accepted and held in abeyance until final commitment is received. Please include a note to this effect on your registration form.

To avoid a LATE REGISTRATION CHARGE, be sure your registration is received prior to NOVEMBER 10th. Late registrations and "walk-ins" at the meeting will be assessed a \$10 LATE FEE.

REIMBURSEMENT of registration fees will be made upon receipt of cancellation notification prior to November 10th. NO REFUNDS WILL BE MADE AFTER THIS DATE. Refunds for unattended workshops, sessions or banquet ticket, WILL NOT be made after arrival to the meeting. Refund will not be made when changing workshop attendance after arrival to the meeting.

For clarification or assistance, please call Roberta Mossdale, NSH Office, 301/552-9678. MAIL registration

and check to: NSH, P.O. BOX 36, LANHAM, MARYLAND 20706.

### NSH/Thomas Edison Program Schedule

**COURSE REVIEWS:** Students interested in review sessions should enroll in one of the workshops scheduled Monday. Workshop registration is free to students formally enrolled in Thomas Edison College. STUDENT MUST SEND PROOF OF ENROLLMENT WITH REGISTRATION FORM.

Tuesday review sessions are only for those planning to take examinations during the convention week, and will allow informal discussion with faculty members, and an opportunity for the examinee to identify and review weak areas. Examinees are encouraged to attend both sessions.

### INTRODUCTORY HISTOTECHNOLOGY/HISTOCHEMISTRY

(Richard Schroeder) HT (ASCP)

Monday - see Workshop #3

Tuesday review, Nov. 17: Powell Room 9 AM - Noon

### HUMAN MICROSCOPIC ANATOMY

(Freida Carson, Ph.D.)

Monday - see Workshop #4

Tuesday review, Nov. 17: Hayden Room 9 AM - Noon

### CURRENT CONCEPTS IN DIAGNOSTIC HISTOPATHOLOGY

(Jules Elias, M.A.)

Monday - see Workshop #5

Tuesday review, Nov. 17: Escalante Room 9 AM - Noon

### EXAMINATION SCHEDULE:

Wednesday, Thursday & Friday, 7-9 AM, BRIDGER ROOM, 2nd floor



**LITTLE  
AMERICA**

"TRADITION  
OF EXCELLENCE"

500 South Main  
Salt Lake City, Utah 84101  
(801) 363-6781

Reservations must be received 2 weeks prior to the opening date of the meeting. Reservations received after 2 weeks cut off will be on a space available basis only. First night's deposit or payment guaranteed by credit card required.

If a room at the rate requested is unavailable, one at the nearest available rate will be reserved.

All rooms subject to applicable city tax.

## NATIONAL SOCIETY OF HISTOTECHNOLOGY

November 16 - 20, 1981

## ROOM RESERVATION REQUEST

NAME \_\_\_\_\_

ADDRESS \_\_\_\_\_

CITY \_\_\_\_\_ STATE/ZIP \_\_\_\_\_

GROUP NAME \_\_\_\_\_

I will be arriving \_\_\_\_\_ AM \_\_\_\_\_ PM on \_\_\_\_\_ DAY \_\_\_\_\_ DATE \_\_\_\_\_

I will depart \_\_\_\_\_ AM \_\_\_\_\_ PM on \_\_\_\_\_ DAY \_\_\_\_\_ DATE \_\_\_\_\_

Credit Card Guar. # \_\_\_\_\_ Exp. Date \_\_\_\_\_ Type of Card \_\_\_\_\_

### ACCOMMODATIONS

Please reserve the following:		Tower
SINGLE OCCUPANCY	\$50	\$65
DOUBLE OCCUPANCY (2 persons, 1 bed)	\$56	\$71
TWIN OCCUPANCY (2 persons, 2 beds)	\$56	\$71

MAIL TO: Little America, 500 South Main, Salt Lake City, Utah 84101



<p>Name _____ (last) _____ (first) _____ (initial)</p> <p>Home Address: _____ (street) _____ (city) _____ (state) _____ (zip)</p> <p>Employer: _____ Address: _____ (street) _____ (city) _____ (state) _____ (zip)</p> <p>Work Telephone No. (_____) _____ (area code) _____</p> <p>Prefer Home _____ Work _____ address for mail.</p> <p>This is your first attendance to an NSH Symposium/Convention. Yes _____</p> <p>Are you an NSH Member? Yes _____ No _____</p> <p>Please check functions you desire to attend and complete Registration.</p>	<p><b>DO NOT USE THIS SPACE</b></p>
---	-------------------------------------

		Workshop Registration	
		Member	Non-Member
		Fees	Fees
#1		\$40	\$45
#2		\$40	\$45
#3		\$40	\$45
#4		\$40	\$45
#5		\$40	\$45
#6		\$40	\$45
#7		\$20	\$25
#8		\$20	\$25
#9		\$20	\$25
#10		\$20	\$25
#11		\$20	\$25
#12		\$20	\$25
#13		\$40	\$45
#14		\$40	\$45
#15		\$40	\$45
#16		\$40	\$45
#17		\$20	\$25
#18		\$20	\$25
#19		\$20	\$25
#20		\$20	\$25
#21		\$20	\$25
#22		\$20	\$25
#23		\$20	\$25
#24		\$20	\$25
Scientific Sessions		\$50	\$60
Banquet		\$18	\$18
<b>\$10 LATE FEE (after Nov. 10th)</b>			
Total Fee: \$ _____			

Workshop Program Monday, Nov. 16			Workshop Program Tuesday, Nov. 17			Wednesday, Thursday & Friday	
	Morning	Afternoon		Morning	Afternoon		
All Day	#7	#10	All Day	#17	#21	Scientific Sessions	
#1	#8	#11	#13	#18	#22	Thursday Banquet	
#2	#9	repeat PM	#14	#19	#23	Wednesday, Open Seminar on Parliamentary Procedure, 7-9 PM; no charge	
#3			#15	#20	#24		
#4	8-11 PM	#12	#16	#20 repeat	PM		
#5							
#6							

## Workshops

### Monday, November 16, 1981

- No. 1: Introduction to Stain Mechanisms**  
*Jack Weiger, HT(ASCP)* 8:30 AM - 4:30 PM  
 Workshop is presented in lecture form with the following subjects discussed in depth: Introduction to electrolytic, self-protein formation; use electric points and how to use them; mordant dyeing; practical applications and various facets of silver reactions, including mechanism of argyrophil and argentaaffin reactions, principles of bacterial staining and mineral reactions. Advanced level workshop.
- No. 2: The Use of Histochemistry in Diagnostic Pathology: An Important Role of Histochemistry in Patient Care**  
*Hugh A. McAllister, M.D., COL, USA* 8:30 AM - 4:30 PM  
 The primary objective of this workshop is to discuss the utilization of special histochemical stains in establishing diagnoses. The role of the histotechnologist in patient care will be illustrated with actual case presentations, as appropriate. Emphasis will be placed on the practical approach to diagnosis of various disease categories.
- No. 3: Introductory Histochemistry/Histochemistry**  
*Richard Schroeder, M.A.* 8:30 AM - 4:30 PM  
 Workshop is designed as an introduction and refresher for the discipline of histochemistry. Concepts of fixation, tissue preparation, sectioning and staining will be presented. Staining procedures utilized as routine in histopathology laboratories, i.e., carbohydrates, lipids, proteins, minerals, bacteria, are discussed. The histochemistry program provides the participant with a more in-depth understanding of routine and sophisticated procedure mechanisms. Cryostat and cryogenic techniques are discussed.
- No. 4: Human Microscopic Anatomy**  
*Freda Carson, Ph.D.* 8:30 AM - 4:30 PM  
 Workshop will emphasize basic tissues and their organization into the different organs. Visual identification and function will be included. Participants will be given a post-test, and handout material will be in a work book format to be filled in during the workshop. Workshop is recommended for Thomas Edison students.

- No. 5: Current Concepts in Diagnostic Histopathology**  
*David Ellis, M.A.* 8:30 AM - 4:30 PM  
 A unique course which encompasses the recommended procedures for handling surgical specimens in the histopathology laboratory. The intimate details of resection as they apply to the handling of small biopsy specimens as well as the rules for grossing tissue from the major human organ systems are included. The course is geared to the experienced histotechnologist working in a clinical environment. It is advisable that the participant purchase a copy of "Basic Techniques in Diagnostic Histopathology" by A. Kennedy, Churchill Livingstone, New York, 1977.
- No. 6: Immunoperoxidase Techniques Utilized in Both Research and Clinical Settings**  
*David Miller, HT(ASCP) & Kay Jenkins, HTL(ASCP)* 8:30 AM - 4:30 PM  
 Limit: 40  
 Workshop will include demonstrations of herpes, toxoplasmosis, viruses, etc. Monoclonal antibodies: Both indirect method and direct methods that are used in a research and clinical setting. Detailed information will be provided in lecture form and handouts.
- No. 7: Principles of Enzyme Histochemistry**  
*Richard W. Denson, Ph.D.* 8:30 AM - 12 Noon  
 Workshop will present an organized treatment of a seemingly bewildering array of procedures. Proper fixation is critical to immobilizing the enzyme in a functioning state. Following that, there are essentially five ways in which a colored product can be produced to mark the site of endogenous enzymes. Direct synthesis, simultaneous coupling, post-incubation coupling, formation and demonstration of azelaic salts, and reduction reactions. Emphasis will be on theoretical aspects so that participants can understand and correct problems experienced with these tests in the laboratory.
- No. 8: Time Management**  
*Renee Nichols, B.A.* 8:30 AM - 12 Noon  
 Limit: 30  
 Developing effective time management is a very important skill that allows you to extend yourself and stretch your abilities. In today's laboratory, time management is a vital skill. At this workshop you will develop forceful plans that will help you learn how to work smarter, not harder, and increase your productivity.







## Workshops, continued

not the true picture, but instead there exists a continuum of slightly different lymphocytes. The latter population of lymphocytes makes a significant contribution to a lymphocyte pool containing at one end of its spectrum T lymphocytes whose surface markers and functions are further removed from those of the B lymphocytes occupying the opposite extreme. Careful consideration of methodological pitfalls and the immunology of T and B cells should be incorporated in any interpretation of enumeration results of these cells in both health and disease.

### 7. Relationships Between Osteoid Widths and Types of Osteoblasts in Bone (A. B. Villanueva, M.A.)

Four types of osteoblasts are normally present lining the intima of osteoid seams. They assume a variety of shapes which we have classified as types 1, 2, 3 and 4 cells. Type 1 cells resemble both mesothelial and fibroblast-like cells frequently present in a very recently

formed osteoid. The osteoid termed here are mostly protruding, collagenic fibers. The nucleus is large, oval shaped and averages 9 micrometers in diameter and 15 micrometers in length, and contains numerous chromatin particles. Tetracycline staining osteoid seams adjacent to these cells are not usually evident. Type 2 cells are the classic, or typical and often designated as active osteoblasts. It assumes a variety of shapes, such as columnar or cuboidal, or pyriform. These osteoblasts can be identified easily by the presence of adjacent nuclear clear zone otherwise known as juxta-nuclear vacuoles or cytoconstrictions, a large, round, clear area at the center of the cell about the size of the nucleus (7um diameter and 10um length) which it adjoins. The width of the osteoid seam adjacent to these cells is approximately 17.4um. Type 3 cells are flattened cells and average about 15 to 20 micrometers in length. The nucleus is spherical or rectangular, centrally or distally located, and averages approximately 4 micrometers in diameter, 10 micrometers in length with an elongated, occasional wavy cytoplasm. The width of the osteoid seam adjacent to these cells is approximately 15.7um. Type 4 cells are transitional, stretched, rectangular-like cells with small, elongated nuclei measuring approximately 2.5 micrometers in diameter and 10 micrometers in length and an extremely thin cytoplasm. The width of the osteoid adjacent to these cells is approximately 5.0um.

## Scientific Sessions

Wednesday, November 18, 1981

### A.M. Session

Malignant Lymphoma, Recent Advances

Open Lung Biopsy

Forensic Histopathology

Lung Cancer — Histology Correlation with Response to Treatment

### P.M. Session

Mortality Study of Histologic Technicians Certified Between 1948-1970

Routine and Special Stains for the Diagnosis of Liver Disease

Estrogen Receptor Analysis: An Immunofluorescent Method

Hormonal Evaluation of the Cytologic Sample

### 7 - 9 P.M.: "Simplifying/Enjoying Parliamentary Procedures"

Session will have a brief introduction to the basics of parliamentary procedure. There will be mock sessions using different situations that illustrate how to accomplish varying results starting at the same place. Finally there will be parliamentary "gamesmanship."

Carl Kjeldsberg, M.D.

Jerrold Abraham, M.D.

Wallace Graham, M.D.

Harmon Eyre, M.D.

Dan Grauman

Randy Lee, M.D.

M. Elizabeth Hammond, M.D.

Catherine Keebler, CT (ASCP) CFIAC

Lois Cook, B.S., HT (ASCP)

Thursday, November 19, 1981

### A.M. Session

Testing for Evidence of Rape

Histochemistry and Sports Medicine

The Pathology of Tumors

### P.M. Session

Alcohol — The Chemical Drug We Drink

Bone Marrow Preparations and Interpretations

Toxicology — Drug Testing

Immunoperoxidase (PAP) Techniques

T. Paulette Sutton, MT (ASCP)

Peter Senzig, HT (ASCP)

Robert Flinner, M.D.

Richard Schroeder, M.A.

Catherine Brunst, HTL (ASCP)

Marilyn Irwin, HT (ASCP) and

Kathy Davis, HTL (ASCP)

Thomas Jennison, Ph.D.

Joe Marty, M.S.

Friday, November 20, 1981

### A.M. Session:

The Interpretation and Application of Diagnostic Special Stains for Surgical Pathology

Animals with a Total Artificial Heart: Blood and Tissue Interactions

Toxic Shock Syndrome

Histologic Technique: Skin Sectioning for the Dermatopathologist

Lecture from Paper Written During Workshop on Tuesday

Erwin Haas, HT (ASCP)

Donald Olsen, D.V.M.

Allen Paris, M.D.

Diane Miller, HT (ASCP)



3. Quickly immerse the cork, with tissue down, into the liquid nitrogen.
4. The frozen specimen is then left in the liquid nitrogen for at least one minute.
5. The cork with frozen specimen attached is then transferred to the cryostat.
6. Wait for a period of 15 to 30 seconds before attaching the cork to a cryostat chuck with OCT embedding compound.
7. Allow the temperature of the frozen specimen to warm to the temperature of the cryostat. Sectioning can now begin.

#### Reference:

1. Moline, S.W., and Glenner, G.G.: Ultrarapid Tissue Freezing in Liquid Nitrogen. *J. Histochem. Cytochem.*, 12:777-783, 1964.

\*OCT embedding compound is a product of LAB-TEK DIVISION, Miles Laboratories, Inc., 30W475 North Aurora Road, Naperville, IL 60566.

†Pre-cut cork discs are available from Slee International, Inc., New York, NY 10011.

## Replies to Inquiry

**Editor's Note:** The following replies were received in regard to why plastic conical tip centrifuge tubes produce more compact buttons than glass conical tip tubes. The original article, authored by Brenda Cuevas, appeared in *Histo-Logic*, Vol. X, No. 2, pg. 146, April 1980.

Priscilla Ann Gregory  
Osteopathic General Hospital  
1750 N.E. 167th Street  
North Miami Beach, Florida 33162

In response to the article by Brenda Cuevas, following is our procedure for Cell Block Preparation.

#### Cell Block Preparation: Plasma - Thrombin Technique

1. Centrifuge the fluid sample at 1500 to 2000 rpm for 15 minutes.
2. Decant supernatant.
3. Add 3 to 5 drops of plasma (plasma with normal PT time).
4. Gently stir with wooden applicator stick to permit plasma to permeate the sediment.
5. Add 3-5 drops of thrombin (Ortho Brain Thromboplastin).
6. Gently stir with wooden applicator stick to allow thrombin to form a clot.
7. After clot forms, place in lens paper and place in embedding cassette.
8. Place in formalin for tissue processing.

#### Reference:

- Pathology Annual, Part I, Vol. 12, University of Miami Medical Center, 1977.

George W. Chang  
University of California  
College of Natural Resources  
Berkeley, California 94720

I was interested to read Ms. Brenda Cuevas' note about "Cell Blocks from Specimens of Body Fluids." I suspect her observation has something to do with the fact that water doesn't really wet the plastic tubes. Thus, there is less water remaining at the bottom of the tube after she pours off the supernatant. The smaller amount of water (fluid) remaining may be less likely to resuspend the compact button of cells.

## Section Adhesive for Paraffin Sections

Alice Esposito  
Methodist Hospital  
Philadelphia, Pennsylvania 19148

The following procedure for adhesion of tissue sections to microscopic slides has been found to be very beneficial in our hands. Some of the benefits are: (1) reduced tissue section detachment from glass slide; (2) improved staining qualities; (3) reduced residue on microscopic slide; (4) less wrinkles evident on finished slide; and (5) saved time.

#### Procedure:

The tissue section ribbon which is obtained from microtomy is placed in the flotation bath containing the following solution. The water bath temperature should be at a minimum of 55°C.

#### Adhesive Solution

Alcohol (ethyl or isopropyl) .....	100.0 ml
Distilled water .....	2400.0 ml
Elmer's Glue-all* .....	6.0 drops
Mix well.	

Sections are picked up on glass slides and placed in a staining rack which is then placed in a slide dryer at 65°C for 15 minutes. Sections are now ready for deparaffinization and staining. Note: Flotation bath and slide dryer temperatures are most important when using this section adhesive. The adhesive solution should be added to water bath just prior to use.

\*Can be purchased at most drugstores and food stores.

## Problems Related to Staining of Connective Tissue Embedded in Water-Soluble Plastics

Peter O. Gerrits  
Department of Anatomy  
University of Groningen  
9713 EZ Groningen  
The Netherlands

In our laboratory we routinely prepare embedding media on the basis of water-soluble methacrylates (JB<sub>4</sub>, GMA, according to Ruddell,<sup>1</sup> 1967, and Sims,<sup>2</sup> 1974). The following problem is experienced.

Although most staining procedures we use in routine histology do not offer any significant problems (occasionally only slight modifications appear to be necessary), we are repeatedly confronted with unsatisfactory results regarding the staining of connective tissue components. Most investigators do not mention similar difficulties, while to our knowledge, an effective connective tissue staining procedure for plastic-embedded specimens has not been reported.

Perhaps someone can provide us with technical protocols or suggestions regarding this issue, in order to obtain sharp contrasts between nuclei and cytoplasm in relation to other tissues. I would gratefully appreciate any suggestions on this matter.

**Editors Note:** Please forward a copy of replies to the Editor, and to Mr. Gerrits at the address above.

#### References:

1. Ruddell, C.L.: Embedding Media for 1-2 Micron Sectioning. *Hydroxethyl Methacrylate Combines with 2-Butoxyethanol*. *Stain Tech.*, 42: 253-255, 1967.
2. Sims B: A Simple Method of Preparing 1-2 Micron Sections of Large Tissue Blocks Using Glycol Methacrylate. *J. Microsc. (Oxf.)*, 101:223-227, 1974.



## MOVED? NEW NAME?

Please print corrections below and return with old label to Lab-Tek. Mailing list updated quarterly. Allow 90 days for address corrections.

Name \_\_\_\_\_

New Address \_\_\_\_\_

City/State \_\_\_\_\_ Zip \_\_\_\_\_



LAB-TEK DIVISION  
Miles Laboratories, Inc.  
30W475 North Aurora Rd.  
Naperville, IL 60540  
Address Correction Requested

BULK RATE  
U.S. POSTAGE  
PAID  
PERMIT NO. 4954  
Chicago, Illinois

**Write on  
and on and on . . . Tissue-Tek III Cassette System™ for accurate,  
time-saving specimen identification.**



The histology systems source

Tired of transferring specimen I.D. numbers with capsules and rings? Now you can have one-time I.D. from grossing to filing, plus three writing surfaces for your permanent specimen identification.

These cassettes have a high resistance to solvents and decalcifying solutions; and they are compatible with microwave techniques. Unique snap-latch and hinge lock provide positive specimen retention and a disposable lid means no more metal covers to clean. All designed to make your lab life easier . . . from LAB-TEK.

Call your LAB-TEK distributor today about the convenient Tissue Tek III Cassette System which includes:

- Uni-Cassette™
- Uni-Cassette in Four Colors for better organization
- Biopsy Uni-Cassette™ with 1mm pores. Helps eliminate the need to wrap most biopsy-type specimens
- Mega-Cassette™ for larger specimens such as bone, eyes and teeth. Special Mega-Cassette base mold also available.

## Can You Help?

Rob Bosma  
Laboratorium voor de Volksgezondheid  
Jelsumerstraat 6, afd. P.A.  
8917 EN Leeuwarden  
The Netherlands

I have a problem in determining and coloring spirochaetes. The methods used are Levaditi's method and the Warthin-Starry method. Who can help me to positively identify (1) *Treponema Pallidum* (syphilis) and (2) *Leptospira Ictrohaemorrhagiae* (or Weil's disease) in liver tissue or lympho-gland tissue (fixed in formalin). Should there be other techniques (i.e., fluorescence methods) for determining of spirochaetes, I would like to be informed.

**Editor's Note:** Please forward a copy of any replies to the Editor, and to Mr. Rob Bosma at the address above.

## An Atlas of Tissue Artifacts

This is the only comprehensive book on tissue artifacts available to histopathology technicians and pathologists. It contains more than 400 photographs illustrating many common and uncommon artifacts produced in all facets of histotechnology. More important, it contains detailed written information on how these artifacts are *identified, produced, prevented and/or eliminated*. This book is the result of more than 20 years of studying tissue artifacts. It is authored by Samuel Wesley Thompson, D.V.M., M.S., and Lee G. Luna, HT (ASCP). The book is available from: Charles C. Thomas, 301-327 E. Lawrence Avenue, Springfield, Illinois 62717.

To receive your own personal copy of HISTO-LOGIC, or to have an associate added to the mailing list, submit home address to: Lab-Tek Division, Miles Laboratories, Inc., 30W475 North Aurora Rd., Naperville, Illinois 60540.  
Printed in U.S.A.

The editor wishes to solicit information, questions, and articles relating to histotechnology. Submit these to: Lee G. Luna, Editor, Histo-Logic, P.O. Box 36, Lanham, Maryland 20706. Articles, photographs, etc., will not be returned unless requested in writing when they are submitted.