

Chloral Hydrate as a Reducer in Silver Impregnation Methods

for the Demonstration of Axons and Neurofibrils

J. R. Goodwin Mater Misericordiae Public Hospital South Brisbane, Queensland, Australia

Bielschowsky type silver impregnation methods are difficult to perform, and there is some difficulty in consistently obtaining good or even adequate results. Bielschowsky's axons and neurofibrils method is cited as being one of the most reliable of the ammoniacal silver impregnation methods. Nonetheless some published procedures carry the following note: "It is common practice with this and other metallic impregnation methods to take several sections from each block. One of these is likely to show rather superior impregnation."

Marsland, et al. developed a method for demonstration of axons in formalin fixed paraffin sections of the brain and spinal cord.¹ This method uses excessive silver nitrate (30 mls of 20% silver nitrate) to prepare the ammoniacal silver solution, and reports of poor staining results as well as difficulty in performing the technique are not uncommon.

A search reveals that all other methods follow the same general path; *i.e.*, an initial silver nitrate impregnation, followed by a reducer (developer). Variations do exist; *i.e.*, some techniques require the impregnation and development steps to be repeated in order to reinforce impregnation of fibers and/or cell processes.

Most silver impregnation methods employ formalin or hydroquinone as a reducer. Attempts to find a more reproducible procedure to demonstrate axons and neurofibrils resulted in the following method, presented here in the hope that it will be as useful to your laboratory as it has been to ours. Microtomy: Paraffin sections

Solutions

10% Silver Nitrate Solution

Silver nitrate	÷.,		÷	÷	÷			ł,	ł.	÷	ł.	÷				ġ,		+	10	g	
Distilled water .			1		-	ļ	2	2	2		į.	1	1	ŝ	1	1	4		100	П	il

3% Sodium Hydroxide Solution

Sodium	hydroxide	****************	3	g
Distilled	water	**************	100	ml

continued on page 236

INSIDE THIS ISSUE

Chloral Hydrate as a Reducer in Silver Impregnation Methods	231
Rapid Dieterle Stain with Microwave Heating	232
A Method for the Removal of Adhered Precipitates from Sections	233
The Immunohistology Spotlight Histo-Tek® Antisera to Bioactive Peptides & Hormones Histo-Tek® Antisera to Cytoskeletal Proteins	234
Immunohistology: Applications in Pathology Program	237

Fixation: 10% buffered neutral formalin

No reader should utilize materials and/or undertake procedures discussed in HISTO-LOGIC articles unders the reader, by reason of education, training and experience, has a complete understanding of the chemical and physical properties of all materials to be utilized and the function of each material by which are procedure is accompliabed copyrelight cases, MILES SCIENTIFIC, DIVESION OF MILES LABORATORIES, INC., NAPERVILLE, ILLINDIS 500545 231

Rapid Dieterle Stain with Microwave Heating

Loretta Sohn and Stanley Shapiro Queens Hospital, Long Island Jewish-Hillside Medical Center Jamaica, New York

Brinn¹ has successfully utilized microwave heating to drastically reduce the staining times of several metallurgical staining techniques; i.e., Grocott's methenamine silver nitrate, periodic acid methenamine silver, Pascual's modification of the Grimelius argyrophil technique, and the Fontana-Masson argentaffin method. A rapid microwave modification of the Dieterle method was in progress, but not sufficiently reproducible at that time for detailed presentation. The application of heat at 60 °C shortened Dieterle's procedure from over 5 hours to 30 minutes.2

Dieterle stain is useful in demonstrating spirochetes in syphilis and Lyme arthritis, bacteria in Legionella pneumophilia and in cat scratch lymphadenopathy, Donovan bodies, and melanin. To demonstrate these entities in the past, our laboratory successfully applied a lengthy modification of the Dieterle method.

Utilizing a modification of the Dieterle technique, accelerated by microwave heating, we have reproducibly demonstrated spirochetes in paraffin tissue sections using only 7 minutes of staining time.

The spirochete stain was free of nonspecific silver deposition and compared favorably with the original long procedure, as well as the Warthin-Starry technique.

Fixation: 10% neutral buffered formalin

Microtomy: Section tissue at 5-6 microns.

Equipment

1. Microwave oven. This laboratory uses a Samsung RE 515 D which brings water to a boil in 65 seconds.

2. Plastic Coplin jar with cover.

Solutions

5% Uranyl Nitrate Solution

Uranyl nitrate	5	g
Ethyl alcohol, absolute	100	ml

10% Gum Mastic Solution

Gum mastic 10 g Ethyl alcohol, absolute 100 ml Allow gum mastic to dissolve in the alcohol for 2 days and filter before use. 1% Silver Nitrate Solution Silver nitrate 1 g 100 ml Distilled water

0.25% Sodium Sulfite Solution Sodium sulfite 0.25 g Distilled water 100 ml Developing Solution

Developing Solution		
Hydroquinone	1.5	g
Sodium sulfite (0.25%)	60	ml
Distilled water	60	ml
Acetone	10	ml
Formaldehyde, concentrated, 37-40%	10	ml
Pyridine	10	ml

Note: Mix ingredients in above order, gently swirling after each addition. The addition of gum mastic turns the solution milky yellow. When the solution is standing in a well lit area, it becomes brown-streaked and should be gently swirled as the streaks appear. The developer is ready for use in about six hours. Store solution at 0-4°C in brown, stoppered bottle. Discard when solution turns dark.

Staining Procedure

- 1. Deparaffinize slides and hydrate to distilled water.
- 2. Place slides in loosely covered plastic Coplin jar containing 5% alcoholic uranyl nitrate solution and heat in microwave oven for 50 seconds.
- 3. Immediately place slide in distilled water at room temperature for 30 seconds.
- 4. Rinse slides in 95% alcohol.
- 5. Place slides in 10% alcoholic gum mastic solution for 3 minutes. The addition of gum mastic turns the solution milky yellow.
- 6. Rinse slides in 95% alcohol.
- 7. Rinse slides in distilled water for 1 minute. Shake slides rapidly and wipe remaining water adjacent to section with gauze.
- 8. Place slides in loosely covered plastic Coplin jar containing 1% silver nitrate solution and heat in microwave oven for 50 seconds.
- 9. Immediately place slides in distilled water at room temperature until cooled.
- 10. Dip slides for 5-30 seconds in developer until sections are pale yellow to light tan.
- 11. Rinse slides in distilled water.
- 12. 95% Alcohol, 2 dips.
- 13. Acetone, 2 dips.
- 14. Xylene, two changes, 15 seconds each.
- 15. Mount coverglass with synthetic mounting medium.

Results

Spirochetes, some bacteria, melanin and Donovan bodies stain black. Background stains pale tan.

Summary

The Dieterle microwave technique is rapid, simple, reproducible and compares favorably with both the original long procedure as well as the Warthin-Starry technique. Nonspecific background silver deposition is absent or negligible.

References:

- 1. Brinn, N.T. (1983) Rapid Metallic Histological Staining Using the Microwave
- Infinit, K.J. (1985) Rapid Methalic Histological Staining Coung the Introview Oven. J. Mistotechnology 6(3), 125.
 Pickett, J.P. and Roggli, V.L. (1982) Rapid Histological Staining Procedures for Material from Immine-suppressed Patients, Am. J. Med. Tech. 11, 893.
 Van Orden, A.E. and Greez, PW. (1977) Modification of the Dieterie Spirocheses Stain, J. Historechnology A(2), 51.
 Wear, D.J., et al. (1983) Cat Scratch Disease: A Bacterial Infection. Science 272: 1405.

221, 1403.

A Method for the Removal of Adhered Precipitates from Sections

Dr. Yuichi G. Watanabe Sapporo Medical College, Sapporo, Japan

In tissue processing laboratories, we often experience troublesome adhesion of precipitates or other unknown material on sections during the course of staining, as described in Histo-Logic[®], Vol 8, No. 2, April 1978. This short article offers a possible solution for this type of problem.

During a staining process, sections of rat submandibular gland were stained in an old solution of Hansen's hematoxylin. Since the staining solution contained a large amount of precipitates, the sections were contaminated as shown in Figure 1. Interestingly enough, precipitates were seen only on a confined area of the tissue; in this case, on the mucous part of the gland.



Figure 2 - Section was decolorized due to the action of the formic acid.

A considerable discoloration of hematoxylin occurred at this step, therefore sections were washed in tap water for a short period and restained.



Figure 1 - Precipitate deposits on mucous part of the gland.

All attempts to remove these precipitates were ineffective, except for one procedure.

The effective method was carried out by dropping a small amount of formic acid (98-100%) on the contaminated sections and rinsing them immediately with water, after which almost all of the precipitates were removed.



Figure 3 - Section after removal of precipitate and subsequent restaining.

Although this method is quite effective, the vigorous action of formic acid occasionally causes partial detachment of sections from slides, especially if they originally adhered incompletely. At present, I do not know how the method works or the extent of its usefulness for removal of precipitates. Absolute ethanol and other solvents have a similar but considerably weaker effect, but they failed to remove precipitates such as those shown herein.



Histo-Tek® Immunohistochemical Grade Antisera

Miles Scientific has offered quality antisera for over 15 years. Our antisera are now being evaluated for suitability in immunohistochemistry.

A sampling of our Histo-Tek[®] immunohistochemical grade antisera follows.

Histo-Tek® Antisera to Bioactive Peptides and Hormones

Histo-Tek antisera to bioactive peptides and hormones are specifically characterized for use in immunohistochemical techniques. Each product is accompanied by a procedural insert with recommended working dilutions. These antisera are tested versus suitable human tissue and tissues from a number of other species and exhibit broad interspecies reactivity.

All antigens are synthetic preparations unless a specific animal source is noted. Antisera to hormones of porcine and sheep origin show strong reactivity to the corresponding human hormone in immunohistological methods.

	Typical	Titers*
	PAP	Indirect FITC
ACTH (1-24)-H	1:2000	1:20
Bombesin-H	1:500	1:20
Calcitonin-H	1:2000	1:20
β Endorphin-H	1:2000	1:20
Gastrin I-H	1:500	1:20
Glucagon-H	1:500	1:20
Growth Hormone (Sheep)	1:2000	1:20
a MSH-H	1:1500	1:20
Neurophysin II (Porcine)	1:5000	1:20
Neurotensin-H	1:500	1:20
Pancreatic Polypeptide-H	1:1000	1:40
Prolactin (Sheep)	1:2000	1:20
Somatostatin-H	1:2000	1:20
Substance P-H	1:1000	1:20
TSH (Sheep)	1:2000	1:20
Vasopressin-Thy	1:2000	1:20

*Incubation – 30 minutes at 37°C. Titers may be lower for shorter or room temperature incubations.

Code No.	Antisera to:	Host Form	Size
64-701-1	ACTH (1-24)-H	rb/ly	0.25 ml
64-702-1	Bombesin-H	rb/ly	0.25 ml
64-703-1	Calcitonin-H	rb/ly	0.25 ml
64-704-1	B Endorphin—H	rb/ly	0.25 ml
64-705-1	Gastrin I-H	rb/ly	0.25 ml
64-706-1	Glucagon-H	rb/ly	0.25 ml
64-707-1	Growth Hormone (Sheep)	rb/ly	0.25 ml
64-708-1	a MSH-H	rb/ly	0.25 ml
64-709-1	Neurophysin II (Porcine)	rb/ly	0.25 ml
64-710-1	Neurotensin-H	rb/ly	0.25 ml
64-711-1	Pancreatic Polypeptide-H	rb/ly	0.25 ml
64-712-1	Prolactin (Sheep)	rb/ly	0.25 ml
64-714-1	Somatostatin-H	rb/ly	0.25 ml
64-715-1	Substance P-H	rb/ly	0.25 ml
64-716-1	TSH (Sheep)	rb/ly	0.25 ml
64-717-1	Vasopressin-Thy	rb/ly	0.25 ml

Thy - Thyroglobulin

H - Hemocyanin

Nomenclatur	e
We have adopted several forms and a all antisers to immunoglobulins are hea pecified to have (H+L) heavy and lig Anti-Mouse IgM should be understoo pecific).	vy chain specific unless ht chain activity, (e.g.,
The following are abbreviations commo tages:	nly appearing on these
ab - antihody ag - antigen Affin Affinity Purified Agarose - Agarose Bound Aggl Agglutinating AP - Alkaline Phosphatase Conj - Conjugate Ferr Ferritin FTTC - Fluorescein FTTC - Fluorescein Prec IgG Fraction Monoform [®] - Monocional HRP - Peroxidase RITC - Rhodamine PBS - phosphate buffered saline SDS - sodium dodecyl sulfate PAGE - polyacrylamide gel electropi	lq - liquid ly - lyophilized dg - dog dk - donkey gp - gainen pig ft - goat mo - mouse rb - mbbit rt - mi sh - sheep



Histo-Tek[®] Antisera to Cytoskeletal Proteins

Miles' Histo-Tek cytoskeletal antisera are polyclonal preparations using highly purified proteins as immunogens. Polyclonal antisera to cytoskeletal proteins show great cross-reactivity with the corresponding cytoskeletal proteins of other mammalian and many non-mammalian species.

They are tested in indirect immunofluorescence on cell cultures and frozen tissue sections from several mammalian and avian species.

All antisera show strong reactivity with human tissues or cells.

	Typical 7	liters*
	PAP	Indirect FITC
Actin	1:50	1:10
Desmin	1:50	1:10
Keratin	1:50	1:10
Myosin	1:50	1:10
Tubulin	1:50	1:10
Also:		
Human Chorionic Gonadotropin (α+β)	1:200	ND
Human Fibronectin	1:200-1:500	1:20
Human Myoglobin	1:200	ND
Human Prostatic Acid Phosphatase	1:200-1:500	1:20-1:50
Porcine Insulin	1:200	ND
*Incubation - 30 minutes at room	n temperature.	

ND - Not Done

Code No.	Antisera to:	Host Form	Size
65-096-1	Actin	rb/lq	1 mi
65-793-1	Deamin	rb/lq	1 ml
65-792-A	Keratin	gp/lq	1 ml
65-791-A	Myosin	rb/lq	1 m
65-095-1	Tubulin	rb/lq	1 m
Also Ava	ilable:	1-1-1	
65-073-1	Human Chorionic Gonadotropin (α+β)	rb/ly	1 ml
64-078-1	Human Fibronectin**	sh/ly	1 m
65-075-1	Human Myoglobin	rb/ly	2 m
65-304-1	Human Prostatic Acid Phosphatase	rb/lq	1 m
65-104-1	Porcine Insulin**	gp/ly	1 m

** - Broad species cross-reactivity.

Ordering Information

Current pricing and terms of sale can be found in the numerical price list published by the Miles Scientific Sales Office serving your country.

Inquiries about product availability or technical data should be made by telephone or letter to your nearest Miles Scientific Sales Office. Ask for Miles Scientific Biochemicals/Immunochemicals.

In countries where it is required, an import license must be secured by the customer prior to placing an order.

Discounts: Standing orders and bulk or multiple package quantities may be eligible for discounts. Contact your Miles Scientific Sales Office for a quotation.

U.S.A.

Miles Scientific

Division of Miles Laboratories, Inc. 2000 North Aurora Road Naperville, IL 60566 Telephone: (800) 348-5465 in Illinois, call (312) 357-3720 aher 5-1-85, (312) 983-5700 Teles: 20037

UNITED KINGDOM EIRE/SCANDINAVIA

Miles Scientific Miles Laboratories, Ltd. Stoke Poges, Slough SL2 4UY Buckinghamshire, England

Telex: 848337 MILABS G WEST GERMANY

Miles Scientific

Telephone: (2814) 5151

Bayer Diagnostik, GmbH Steinerstrasse 15 PO. Box 701 360 D-8000 Munich 70 Telephone: 089-724030 Teles: 522406 MILES D

AUSTRALIA NEW ZEALAND

Miles Scientific Miles Laboratories Australia (Pty), Ltd. Private Bag 5 Mulgrave North 3170, Victoria Telephone: (03) 560-0611 Teles: MILES AA31129

JAPAN

Miles Scientific Miles Labs. Inc.-Japan Br. Aoyama Bidg. 2-3. Kha-Aoyama 1-chome Minato-ku, Tokyo 107 Telephone: (03) 403-9951 Teles: 242-3548

CANADA

Miles Scientific

Miles Laboratories, Ltd. 77 Belfield Road Rexdale, Ontario M9W 1G6 Telephone: (416) 248-0771 Teles: 06-989399

ITALY/SWITZERLAND GREECE

Miles Scientific

Miles Italiana, S.p.A. via F.L. Miles 10 20040 Caremago Brianza Milano Telephone: 02-95019231 Telex: 332016 MILESCAB 1

FRANCE/BENELUX

Miles Scientific

Laboratories Miles S.A. Tour Maine Montparnasse 33 avenue da Maine 75755 Paria – Codex 15 Felephone: Paris 5385243 Telex: MILAB 204315 F

SPAIN/PORTUGAL

Miles Scientific Miles-Martin Laboraturies, S.A.E. Plaza de Espana, 10

Madrid-13 Telephone: Madrid 242-5104 Teles: 22590 MILAB E

continued from page 231

Ammoniacal Silver Solution

To 5 ml of 10% silver nitrate add just enough ammonium hydroxide, drop-by-drop, to redissolve the precipitate which forms. Add 5 ml of the 3% sodium hydroxide solution and with a few drops of ammonium hydroxide redissolve any silver precipitate which may form.

5% Chloral Hydrate Solution

Chloral hydrate	• •	14	 		 4	 +	2		 4	5 g	1
Distilled water			 	 	 	 		 		100 n	nl

5% Sodium Thiosulfate Solution

Sodium	thiosulfa	te			+	+	÷	÷	é	-		÷	÷		÷	-	4		+	5	g	
Distilled	water .	+	2	-			2	÷	-	-	ŝ			-		ł			1	100	ml	

0.02% Gold Chloride Solution

Gold chloride	0.02	g
Distilled water	100	ml

5% Oxalic Acid Solution

Oxalic acid	 5 g
Distilled water	 100 ml

Aniline Blue Solution

Aniline blue	0.1	g
Oxalic acid	5	g
Distilled water	500	ml

Staining Procedure

- 1. Deparaffinize slides to the last absolute alcohol.
- From this solution, celloidinize slides and harden in 70% alcohol.

Note: To produce even staining across the entire slide, a smooth even coat of celloidin is essential. To do this, dip slide in celloidin, drain to one corner and air dry for 4-5 seconds before immersing in 70% alcohol.

 Place slides in 10% silver nitrate solution at room temperature for 1 hour.

Note: The initial 1 hour in 10% silver nitrate is necessary to successfully stain neurofibrils. If only axons are required, the time can be shortened to 10 minutes. Overnight staining in a less than 10% solution of silver nitrate has not been tried and may result in less background staining. Background staining produced by this method is not excessive, although it would be advantageous to eliminate it.

- 4. Wash slides well in distilled water.
- Place slides in ammoniacal silver solution for 2 minutes.
- 6. Wash slides well in distilled water.
- Place slides in 5% chloral hydrate solution for 1 minute.
- 8. Wash slides well in distilled water.
- Place slides in second change of ammoniacal silver solution for 2 minutes.

 Wash off ammoniacal silver with 5% chloral hydrate and leave 1 minute.

Note: Pour enough 5% chloral hydrate on the slide so that it approximately equals the amount of ammoniacal silver on the slide. This solution turns black. Leave for approximately 15 seconds and completely wash off with chloral hydrate.

 Place slides in 5% sodium thiosulfate solution for 1 minute.

Note: This technique may have an additional role. While the sodium thiosulfate is needed to remove unreduced silver, it also produces a blackening of the celloidin. However, it is difficult to observe whether or not this blackening contributes to the staining of the axons and neurofibrils.

- Wash slides well in distilled water and remove celloidin.
- Place slides in 0.02% gold choride solution for at least 10 minutes.

Optional: If only axon staining is required.

Note: Without the toning step, black staining axons on a light brown to brown background are produced. However, to properly delineate neurofibrils, toning is necessary. Stronger gold chloride solutions have a deleterious effect on the clarity of axon and neurofibril staining.

- 14. Rinse slides with distilled water.
- Optional: 5% oxalic acid for about 15 minutes. This step reddens the background and produces an intensifying effect on the stained fibers.
- Optional: Pour off oxalic acid and apply aniline blue solution for 5 minutes, then quickly rinse in tap water. This step is necessary if further staining of the background is required.
- Dehydrate in alcohol, 95%, absolute, and clear in xylene, three changes each.
- Mount coverglass with resinous media.

Results



Figure 1

continued on page 238



Now we are about to establish a new tradition: Educational support for emerging and developing technologies.

We're announcing the first of our educational programs for 1985:

Immunobistology: Applications in Pathology

This program is designed for both Pathologists and Histotechnologists. The course deals with the practical use and implementation of immunohistology techniques. It will combine seminars by noted scientists in the field, a discussion of DRG's and other related economic and regulatory issues, as well as a case consultation workshop and a hands-on wet workshop.

The course is jointly presented by University MD-PhD Educators and Miles Scientific.

Miles Scientific



2000 North Aurora Road, Naperville, IL 60566

To learn more about this and other educational programs in your area, please fill out the form below and return it to Miles Scientific.

Mail to:	Miles Scientific 2000 North Aurora Road Naperville, IL 60566 Attn: Educational Program Coordi	nator
Please se	end information on educational prop	grams in my area.
Name:		
Title:		
Institutio	on:	
Address:		
City/State: Zi		Zip:
Phone:		

continued from page 236





Figure 3

Figures 1, 2 and 3 are representative of the results possible with this method. Note the extensive demonstration of fine, delicate neurofibrils in figure 1 and 3.

Axons and neurofibrils — black Some connective tissue and elastica — black Background — reddish or grey

This method produces clean, crisp slides of a nerve fiber stain in a reasonable period of time. The method is consistently reproducible with economical use of silver. We trust this method will be welcomed as a replacement for the procedures you now use.

Reference:

 Marsland, T.A., et al. (1954) Modification of the Gless Silver Impregnation for Paraffin Sections, J. Neuropathology and Experimental Neurology 13, 587.

Miles Scientific Division of Miles Laboratories, Inc. 2000 North Aurora Rd. Naperville, IL 60566

Address Correction Requested



To receive your own copy of HISTO-LOGIC, or to have someone added to the mailing list, submit home address to: Miles Scientific, Division of Miles Laboratories, Inc. 2000 No. Aurora Rd., Naperville, IL 60566. Printed in U.S.A. The editor wishes to solicit information, questions, and articles relating to histotechnology. Sabmit these to: Lee G. Luna, Editor, Histo-Logic, PO. Box 36, Lanham, Maryland 20706, Articles, photographa, etc., will not be returned unless requested in writing when they are submitted.