


Tissue-Tek® Archive Systems



SAKURA

Tissue-Tek® Lab Aid®

Storage Cabinets

Safe, Secure, and Trusted...

Tissue-Tek® Lab Aid® Slide Filing System

Trusted for more than 60 years, the Tissue-Tek® LabAid® allows for convenient storage and retrieval of microscopic slides or transparencies. Two different vertical filing drawer sections can be stacked in any combination. 1" drawer section, 14 drawers holds 465 (1" x 3") slides. 2" drawer section, 7 drawers holds 271 (3" x 2") slides or 271 (2" x 2") transparencies.



- Requires only 19 x 19 inches of floor space
- Backstops easily attach inside drawers to hold slides in place
- Welded steel for stackability and durability
- Product Codes: 4010 thru 4027

Dimensions: 1" x 14" drawer or 2" x 7" drawer

Tissue-Tek® Lab Aid® Ultra

The Tissue-Tek® LabAid® Ultra provides a large storage capacity for slides and paraffin blocks. It consists of ten smooth-rolling drawers for storing over 81,900 slides or 21,000 blocks. Constructed from sturdy steel for stackability and durability.



- Drawer lock prevents tampering and ensures specimen safety
- Slide or block trays can be easily removed and conveniently repositioned
- Foam inserts are available for partially filled trays to secure slides and blocks in the proper upright position
- Slides and blocks can be stored in the same cabinet
- Product Code: *Slide Cabinet: 3991 Block Cabinet: 3992*

Dimensions: 23.6"(w) x 29.5"(d) x 43.3" (h), weight: 275 lbs

Tissue-Tek® Block Filing Cabinet

The Tissue-Tek® Block Filing Cabinet enables organization and storage of blocks.



- One piece construction eliminates hinges and doors for dust free storage
- Six impact resistant plastic drawers feature generous space for identification
- Stackable
- Available with or without Tissue-Tek® Uni-Cassettes®
- Product Code: Cabinet only : 4192

Dimensions: 17"(w) x 9.5"(d) x 15"(h)



Trusted Innovation™

Sakura Finetek USA, Inc.
1750 West 214th Street
Torrance, CA 90501
Phone: 800-725-8723
www.sakuraus.com