HIST@-LOGIC®

No resider should utilize materials and/or undertake procedures discussed in HISTO-LOGIC articles unless the reader, by reason of education, training and experience, has a complete understanding of the chemical and physical properties of all materials to be utilized and the function of each material by which any procedure is accomplished.

Editor, Lee G. Luna, D. Lit., H.T. (ASCP)

Technical Bulletin for Histotechnology Published: January, April, July, October

Vol. VIII, No. 3 - July, 1978

Golden Forceps Award Winner

We are pleased to announce that Anna Marie England has been selected as the recipient of the Golden Forceps Award for 1978. Her paper, "Gram-Positive and Gram-Negative Bacteria Control," was selected from articles submitted to Histo-Logic during the past year. Criteria for



selection are clarity, originality and scientific contribution. The Golden Forceps Award will be presented at the Symposium/Convention of the National Society for Histotechnology to be held in Oklahoma City, Oklahoma, September 25-29, 1978. Reprints of the winning article which appeared in the January, 1978 issue of *Histo-Logic* are available from Lab-Tek Division, Miles Laboratories Inc., 30W475 North Aurora Road, Naperville, Illinois 60540.

Gelatin-Chrome Alum: A Better Section Adhesive

Cathy Adler* St. Vincent Charity Hospital Cleveland, Ohio 44115

Mayer's egg albumen is the traditional section adhesive which has been used extensively for many years. Its value as a section adhesive is well recognized, particularly if one sections brain, keratinized material and decalcified bone. These and other similar materials have a tendency to fall off the slides during the staining process and therefore require an effective adhesive.

According to Lillie, the effectiveness of this solution is probably due to its action as a surface tension depressant, since the constituents are readily soluble in water. Egg albumen residue often stains and this produces an objectionable artifact on completed slides and tissue sections. In addition, slides on which egg albumen has been used must remain in the drying oven somewhat longer to insure effective adhesion.

A comparative study of four section adhesives was conducted to determine which of the four proved the most effective and eliminated the artifacts produced on slides by the use of egg albumen. The test sections were obtained from large brain specimens and decalcified dense bone specimens. Duplicate sections were obtained from each block and mounted on slides which contained the different adhesives. Test slides were dried in a 58° C oven for 45 minutes, 2 hours, and 24 hours.

ADHESIVES USED FOR STUDY Mayer's Egg Albumen

White portion of an egg Glycerin

Equal Parts

Loosen egg white with a fork before mixing with glycerin.

Mix both ingredients well and filter solution through several layers of gauze. Add a crystal of thymol to the solution to prevent bacteria and fungal growth.

Starch Paste²

Powdered starch 1.0 gm
Cold water
Mix well and pour this solution into 30.0 ml of boiling water.
Add 2 drops of 1N hydrochloric acid. Boil solution for 5 min-
utes. Add a crystal of thymol.

5% Gelatin

Gelatin														5.0	gm
Cold water					*			+						15.0	ml

Allow this solution to stand for a short period until the gelatin becomes soft. Add 85.0 ml boiling water to above solution. Continue gentle heat with bunsen burner or hot plate if gelatin is not completely dissolved.

Gelatin-Chrome Alum²

Gelatin		+		+	4		+		+		+	į				1.0	gr	n
Cold water		1					4			·						5.0	m	1

Allow this solution to stand for a short period until the gelatin becomes soft. Add 85.0 ml boiling water to above solution. Continue gentle heat with bunsen burner or hot plate if gelatin is not completely dissolved. Add 0.1 gm of chrome alum to the solution and mix well.

MODE OF ADHESIVE APPLICATION TO SLIDES

Mayer's Egg Albumen: Apply evenly on clean glass slide with a finger or small sponge.

Starch Paste: Same as egg albumen solution.

NOTE: This adhesive should not be used if periodic acid Schiff type procedures are to be performed on sections since the starch will react with this procedure.²

5% Gelatin: Add 10.0 ml to flotation bath and stir well. Gelatin-Chrome Alum: Same as 5% gelatin solution.

Results:

The brain sections remained attached to the slides with all adhesives during the staining process. However, some stain residue (in varying degrees) was observed with all adhesives with the exception of the gelatin-chrome alum. Bone sections which were dried in a 58° C oven for 45 minutes fell off the slide during the staining process with all adhesives except the gelatin-chrome alum. The only bone sections remaining attached to the slides with 5% gelatin were those which were dried in a 58° C oven for 24 hours. All other sections using this adhesive and various drying periods fell off the slides. Bone sections placed on the starch paste adhesive and the Mayer's egg albumen remained on the slide only when exposed for 24 hours drying in the 58° C oven.

Conclusion:

The evaluation of these four section adhesives has convinced us that the gelatin-chrome alum solution is less likely to allow sections to detach from the glass slide during the staining process. In addition, many special stains, including silver reactions, can be performed without any chemical interference or accumulated debris.

*Ms. Adler was a histology student at St. Vincent Charity Hospital when she submitted this article for publication.

References:

- Lillie, R.D.: Histopathologic Technic and Practical Histochemistry, 3rd ed., p. 77, McGraw-Hill, New York, 1965.
- Culling, C.F.A.: Handbook of Histopathological and Histochemical Techniques, 3rd ed., p. 131, Butterworths, Reading, Mass., 1974.

A Modified Thiosulfation-Aldehyde Fuchsin Method for Pancreatic Beta Cell Granules

Jack B. Wenger and Lee G. Luna Washington, D.C. 20306

Introduction:

In 1974 Bussolati and Bassa¹ reported the use of a solution of ammonium hydroxide and cupric sulfate as a pretreatment to the aldehyde fuchsin method for the demonstration of pancreatic islet B cells. They found this treatment, which they referred to as "thiosulphation," provided more specific and darker staining of beta cell granules. We have modified the procedure by substituting ammonium sulfate [(NH₄)₂ SO₄] for ammonium hydroxide and dividing the pretreatment solutions into individual stock solutions. These are stable solutions, while the ammonium hydroxide-cupric sulfate-sodium metabisulfite solution used by Bussolati and Bassa is somewhat unstable.

In addition, we have substituted potassium dichromate for paraldehyde in the aldehyde fuchsin solution. This was done for two reasons:

- 1. Paraldehyde is a narcotic and therefore a controlled item.
- 2. This chemical is often difficult to obtain.

This procedure has been used in our laboratories for over a year with excellent results. It has provided us with darker and more specific demonstration of pancreatic B cell granules.

Fixation: Buffered neutral formalin

Microtomy: Cut paraffin sections at 6 micra.

Solutions:

Thiosulfation Solution A (Stock)

Ammonium sulfate [(NH ₄) ₂ SO ₄]6.0 gm	1
Copper sulfate, hydrate (Cu SO ₄ 5H ₂ O) 2.5 gm	1
Distilled water	

Thiosulfation Solution B (Stock)

Sodium hydroxide.								2.0 gm
Distilled water					4	+		400.0 ml

Thiosulfation Solution C (Stock)

Sodium	meta-bisulfite	m
Distilled	water	nl

Solutions A and B keep well, provided evaporation is minimized. Solution C remains stable for about one week.

Thiosulfation Solution (Working)

Mix in the following order immediately before use:

Solution	Α.		+	4	*		4							-		-		20	.0	m	ı
Solution	B.			ï		U	ı	Ċ	ı			ı	ı		U	ď	3	20	0	m	ı

Solution is now dark blue due to the complexing of the ammonia and copper (pH is approximately 9).

ADD:

Solution	C			5.0 ml
Solution turns	bluish-green	and may	have a	light precipitate
(solution will b	ecome nale an	d clear as	the reacti	on progresses)

Wenger's Aldehyde Fuchsin

Basic fuchsin		+		+	+		į.	·	+		2.5	gm
Ethyl alcohol, 70%*			,	,		,		,	,	,	500.0	ml
Hydrochloric acid				1					÷		5.0	ml
Potassium dichromate												

Combine the above ingredients and cap the container tightly. Allow solution to set at room temperature for five days. Check daily for the formation of acetaldehyde by smelling solution.**

Acetaldehyde can be detected by its characteristic pungent fruity odor. If acetaldehyde is not detected within 24 hours, ethyl

*Ethyl alcohol is necessary due to its reactivity with the potassium dichromate to form the acetaldehyde. alcohol and/or potassium dichromate was omitted, or possibly an insufficient amount of hydrochloric acid was used. At the end of the five day period, check the quality of the solution by staining a known elastic control. If the elastic tissue is well demonstrated, the solution is good and should be placed in a refrigerator.

0.25% Metanil Yellow

Metanil yellow	4.		4				4		+			+		0.25	gm
Distilled water													٠	100.0	ml
Glacial acetic a		-9												0.25	gm

Staining Procedure:

- 1. Decerate and dehydrate slides to distilled water.
- Place slides in fresh working thiosulfation solution for 2 hours.
- Wash slides in tap water for 1 minute and rinse in distilled water.
- 4. Rinse slides in 95% alcohol
- Place slides in Wenger's aldehyde fuchsin solution for 2 minutes.
- 6. Rinse slides in 95% alcohol until no more stain comes off the

CHECK CONTROL SLIDE MICROSCOPICALLY. IF ISLET CELL GRANULES ARE PALE OR UNSTAINED, REPEAT STEPS 2 THROUGH 6 WITH FRESH SOLUTIONS.

- 7. Counterstain slides in metanil yellow for 1 minute.
- Dehydrate in 95% alcohol, absolute alcohol, and clear in xylene, 2 changes each.
- 9. Mount coverglass with synthetic resin.

Results

Beta cell granules - purple Background - tan to yellow

Reference:

 Bussolati, G. & Bassa, T.: Thiosulfation Aldehyde Fuchsin (TAF) Procedure for the Staining of Pancreatic B Cells. Stain Tech., Vol. 49, No. 5, pp. 313-315, 1974.

Patron Program

National Society for Histotechnology

The National Society for Histotechnology has instituted a Patron Program to allow those individuals who are interested in the Society's activities to support its many programs with monetary contributions. These contributions are tax deductible.

The program functions in the following manner. Patrons are asked to make checks payable to the National Society for Histotechnology and indicate on the check, or an attached note, that the contribution is for the Patron Program. NSH will, upon receipt of contribution, forward a receipt to patron(s) acknowledging amount of contribution (for tax purposes). The patron's full name and amount of contribution will be printed in the Journal of Histotechnology once a year.


The contribution can be in any amount, but it is recommended that no contribution be less than \$5. Contributions can be made as often as one wishes: monthly, bimonthly, quarterly, or once a year. It is important that potential patrons understand they are not committed to continual contributions unless they choose to do so. Mail contributions to: National Society for Histotechnology, P.O. Box 36, Lanham, Maryland 20801.

Examination of the Central Nervous System (Spinal Cord) in Situ

M. E. Shelton, D.V.M., Rose M. Aspengren and Grace Faber College of Veterinary Medicine Iowa State University Ames, Iowa 50011

The use of a 10-15% (saturated) solution of disodium EDTA* for decalcification has been well recognized. 1.2.3 An advantage to its use as a decalcifying agent is that the stainability of bone and

^{**}Use care in smelling solution.



COLUBITION

TISSUE-TEK® II PROCESS/EMBEDDING CASSETTES NOW AVAILABLE FROM YOUR LAB-TEK DISTRIBUTOR IN FOUR DISTINCTIVE COLORS:

YELLOW, PINK, GREEN AND BLUE

#4179

#4182

#4183

#4184

Identical in design to 4187 Tissue-Tek II Process/Embedding Cassettes (white), color cassettes feature a slanted writing surface for quick, positive identification. Perforated cassette and process cover insure thorough and rapid penetration of processing fluids. Color Cassettes are compatible with tissue processing solutions used for histological preparations.

VERSATILE COLOR CASSETTES OFFER ENDLESS POSSIBILITIES TO FIT YOUR NEEDS

- Separate surgicals from autopsies
- · Alert staff to stat biopsies
- Identify important specimens for retention and/or teaching purposes
- Identify frozen sections
- Separate specimens which require different processing or special techniques
- · Differentiate satellite hospitals using central facilities
- · Segregate cytology specimens
- Isolate research specimens
- Segregate bone marrow and decalcified specimens.

ANOTHER WAY
TISSUE-TEK
SYSTEMS
SERVE YOU FROM
START TO STORAGE.



LAB-TEK...setting the standards by which performance is judged.

NOW IT'S EASY TO CREATE YOUR OWN COLOR IDENTIFICATION SYSTEM WITH TISSUE-TEK® II PROCESS/EMBEDDING CASSETTES

BUSINESS REPLY MAIL

NO POSTAGE STAMP NECESSARY IF MAILED IN THE UNITED STATES

PERMIT NO. 153 NAPERVILLE.

ILLINOIS 60540

POSTAGE WILL BE PAID BY

LAB-TEK DIVISION

Miles Laboratories, Inc. 30W475 North Aurora Road Naperville, IL 60540

Attn: Marketing Services

TO LEARN MORE

about this new color identification system, simply fill out the coupon below, refold, staple and return to us. No postage needed.

Name/Title______
Institution _____
Address

City/State/Zip

I presently use: ☐ 4151 Tissue-Tek Embedding Ring ☐ 4187 Tissue-Tek II Process/Embedding Cassette (clear)

associated soft tissue is not adversely affected following immersion in the decalcifying solution for several weeks or months. Two disadvantages to its use in histological preparations are predominant. Because of its slowness of action, decalcification may not be complete for several weeks. Material decalcified for too long (over 5 months) may exhibit diluted staining reactions.

In the study of the immature spinal cord many problems are encountered. One problem is the poor fixation, or more specifically, the softness of the tissue after perfusion and fixation in 10% buffered neutral formalin (BNF). The inability to prepare this soft tissue for embedding without gross artifact made a new way for examination imperative. After perfusion through the carotid with 10% BNF, the vertebral column was removed with the spinal cord and dorsal root ganglion in situ. The column was bisected through the twelfth intervertebral disc and placed in fresh 10% BNF. After allowing time for proper fixation, the tissue was cut between vertebral bodies for paraffin embedding, The tissues were placed in individual polyethylene plastic tissue capsules.** (This is important in that the stainless steel capsules seem to deteriorate after much use in this solution.) Complete vertebral columns from adult rodents were immersed and cut for blocking after decalcification. After two to five months in the EDTA solution, the tissue cubes were embedded and sectioned at 6 and 10 microns. The solution was changed every two months. H & E and cresyl echt violet stains were used. Stainability of bone and soft tissues was excellent. Adult rodents, decalcified for 5 months, and young dogs which were less than eight weeks of age, decalcified for two months, have been processed this way (Figures 1-5).




FIGURE 1:

Adult mouse spinal cord at low lumbar level. Whole vertebral column fixed in 10% BNF, decalcified in 10% (saturated) sodium EDTA for 5 months. Animal had spina bifida. (Cresyl echt violet, x 10, 10 microns)




FIGURE 2:

Upper sacral level of same vertebral column as Figure 1. (Cresyl echt violet, x 10, 10 microns)




FIGURE 3:

First lumbar section from 2 day old canine. Animal was perfused with 10% BNF. Decalcification proceeded for 2 months in 10% sodium EDTA. (Cresyl echt violet, x 4, 10 microns)




FIGURE 4: Section from sixth lumbar level from same animal as Figure 3. (H & E, x 4, 10 microns)




FIGURE 5:

Section from upper sacral level. Same as Figure 3. (Cresyl echt violet, x 10, 10 microns)

References

- Wallington, E.A.: Histological Methods for Bone. Butterworth and Co., Reading, Mass., 1972.
- Brain, E.G.: The Preparation of Decalcified Sections. Charles C Thomas, Springfield, Ill., 1966.
- 3. Graham, D.: Personal communication.
- *Ethylenediaminetetra-acetic acid, disodium salt: C₁₀H₁₄O₈N₂Na₂, 2 H₂O
- **TC-3 Plastic tissue capsules, Lipshaw Manufacturing Corp., Detroit, Michigan 48210.

HISL®-FOCIC.

Plan to altered NSH Sympoolium/Convention Sept. 25-29 in Oblahoma City




AeT-de.J Division

Address Correction Requested

30W475 N. Aurora Rd., Naperville, IL 60540

Please print corrections below and return with old label to Lab-Tek. Mailing list updated quarterly. Allow 90 days for address corrections.



Tissue-Tek® II Process/Embedding Cassette is numbered at the start for positive identification. From start to storage, tissue and cassette travel together to eliminate misidentification. And the cassette does double duty: first it serves as a processing capsule, and then as an embedding ring. It forms streamlined tissue blocks that save 30-50% on paraffin costs and take up 50% less space than previous methods. Available in clear plastic No. 4191 or the highly resistant white plastic No. 4187. Process Cover No. 4194. The polished steel base mold makes a clearer casting which presents a superior cutting face.

Its clean sharp edges insure an excellent fit in the microtome. The mold's design also helps reduce cleanup, storage, and housekeeping chores. Base molds available in five sizes. Nos. 4161-4165.

Lab-Tek . . . setting the standards by which performance is judged.

Technically Speaking

PERMIT NO, 4954 Chicago, Illinois

CIIVA

U.S. POSTAGE

BULK RATE

Robert A. Clark, Technical Services Lab-Tek, Naperville, IL 60540

Lab-Tek has received numerous inquiries regarding the installation and use of the exhaust deflector used on Tissue-Tek II Tissue Embedding Center.

Like any true refrigeration system, the compressor will operate more efficiently

if the warm air generated by the compressor is diverted away from the unit.

The exhaust deflector provided with the embedding center has several functions. It directs warm air generated by the compressor away from the unit and prevents the unit from being placed directly against a wall.

If your deflector has been lost or misplaced, it can be ordered from any authorized Lab-Tek distributor. The part number is 50090021.



- 1. Turn unit around so that rear of unit faces you.
- 2. Loosen screw which affixes paraffin reservoir to base unit.
- Examine foil vent for accumulation of lint and dust. If necessary, remove vent, wash in hot, soapy water, rinse, dry and replace.
- 4. Locate round opening in foil vent.
- Position exhaust deflector so that large scoop faces up and is centered over round opening. Make sure that left side of exhaust deflector extends approximately 1/8" over trim panel of base unit.
- Holding exhaust deflector against rear of embedding center, compress both sides so that left lip of deflector slips under trim of base unit. Release pressure.
- Lift rear portion of paraffin reservoir so that slot (located top right) engages under lip of paraffin reservoir.
- Move paraffin reservoir so that it is properly positioned on base unit. Replace screw removed in step #2.

If additional information is needed, please write Lab-Tek Customer Service.

To receive your own personal copy of HISTO-LOGIC, or to have an associate added to the mailing list, submit home address to: Lab-Tek Division, Miles Laboratories, Inc., 30W475 North Aurora Rd., Naperville, Illinois 60540.

Printed in U.S.A.

The editor wishes to solicit information, questions, and articles relating to histotechnology. Submit these to: Lee G. Luna, Editor, Histo-Logic, P.O. Box 36, Lanham, Maryland 20801. Articles, photographs, etc., will not be returned unless requested in writing when they are submitted.